

Product Datasheet

CD117/c-kit Antibody (YB5.B8) - BSA Free NBP1-43358-0.1mg

Unit Size: 0.1 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-43358

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-43358



NBP1-43358-0.1mg

CD117/c-kit Antibody (YB5.B8) - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	5 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	YB5.B8
Preservative	0.09% Sodium Azide
Isotype	IgG1 Kappa
Purity	Protein A or G purified
Buffer	PBS (pH 7.2)


Product Description	
Description	Novus Biologicals Mouse CD117/c-kit Antibody (YB5.B8) - BSA Free (NBP1-43358) is a monoclonal antibody validated for use in IHC, Flow and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	3815
Gene Symbol	KIT
Species	Human
Reactivity Notes	This antibody is reactive to Human.
Marker	Hematopoietic Stem Cell Marker
Immunogen	Human CD117

Product Application Details	
Applications	Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation
Recommended Dilutions	Flow Cytometry 1:10-1:1000, Immunohistochemistry 1:10-1:500, Immunoprecipitation 1:10-1:500, Immunohistochemistry-Frozen 1:10-1:500
Application Notes	The YB5.B8 antibody has been tested by flow cytometric analysis of the TF-1 cell line, a human erythroleukemia cell line. This can be used at less than or equal to 1 ug per test. Cell number should be determined empirically but can range from 10 ⁵ to 10 ⁸ cells/test.




Images

Flow Cytometry: CD117/c-kit Antibody (YB5.B8) [NBP1-43358] - Analysis using the DyLight 405 conjugate of NBP1-43358. Staining of CD117 in human PBMCs using anti-CD117 antibody. Image from verified customer review.



Flow Cytometry: CD117/c-kit Antibody (YB5.B8) [NBP1-43358] - Analysis using the Biotin conjugate of NBP1-43358. Staining of TF-1 cell line with 0.5 ug of Mouse IgG1 kappa Isotype Control Biotin (open histogram) or 0.5 ug of Anti-Human CD117 (c-Kit) Biotin (filled histogram) followed by Streptavidin PE.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-43358-0.1mg

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
NBP1-43358V	CD117/c-kit Antibody (YB5.B8) [DyLight 405]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-43358

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

