

Product Datasheet

CD27/TNFRSF7 Antibody (LG.7F9) - BSA Free NBP1-43428-0.1mg

Unit Size: 0.1 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-43428

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-43428



NBP1-43428-0.1mg

CD27/TNFRSF7 Antibody (LG.7F9) - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	LG.7F9
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Protein A or G purified
Buffer	PBS (pH 7.2)


Product Description	
Description	Novus Biologicals Armenian Hamster CD27/TNFRSF7 Antibody (LG.7F9) - BSA Free (NBP1-43428) is a monoclonal antibody validated for use in Flow and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Armenian Hamster
Gene ID	939
Gene Symbol	CD27
Species	Human, Mouse, Rat
Reactivity Notes	This antibody is reactive to Mouse, Rat, Human.
Immunogen	The immunogen for this antibody was CD27.

Product Application Details	
Applications	Flow Cytometry, Immunoprecipitation
Recommended Dilutions	Flow Cytometry 1:10-1:1000, Immunoprecipitation 1:10-1:500
Application Notes	The LG.7F9 antibody has been tested by flow cytometric analysis of mouse splenocyte cell suspensions. This can be used at less than or equal to 0.5 ug per test. Cell number should be determined empirically but can range from 10 ⁵ to 10 ⁸ cells/test.




Images

Flow Cytometry: CD27/TNFRSF7 Antibody (LG.7F9) [NBP1-43428] - Analysis using the FITC conjugate of NBP1-43428. Staining of C57BL/6 splenocytes with Anti-Human/Mouse CD45R (B220) APC and 0.25 micrograms of Armenian Hamster IgG Isotype Control FITC (left) or 0.25 micrograms of Anti-CD27 FITC (right).



Flow Cytometry: CD27/TNFRSF7 Antibody (LG.7F9) [NBP1-43428] - Staining of mouse splenocytes with 0.25 ug Hamster IgG Isotype Control Purified (open histogram) or 0.25 ug Anti-CD27 Purified (filled histogram) followed by Anti-Armenian Hamster IgG FITC . Total viable cells were used for analysis.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-43428-0.1mg

NB120-5738	Goat anti-Armenian Hamster IgG (H+L) Secondary Antibody
NBP1-97046	Armenian Hamster IgG Isotype Control
NBP1-43428APC	CD27/TNFRSF7 Antibody (LG.7F9) [Allophycocyanin]
NBP2-38434PEP	CD27/TNFRSF7 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-43428

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

