

Product Datasheet

CD3 zeta Antibody (H146-968) [PE] NBP1-44683PE

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-44683PE

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-44683PE



NBP1-44683PE

CD3 zeta Antibody (H146-968) [PE]


Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	H146-968
Preservative	0.05% Sodium Azide
Isotype	IgG2 Kappa
Conjugate	PE
Purity	Protein A purified
Buffer	PBS

Product Description	
Host	Armenian Hamster
Gene ID	919
Gene Symbol	CD247
Species	Human, Mouse
Specificity/Sensitivity	The Armenian hamster antibody H146-968 reacts with CD3 zeta chain (CD247), which is a component of TCR/CD3 complex expressed on T cells.
Immunogen	Synthetic peptide corresponding to amino acids 151-164 of mouse CD3 zeta.

Product Application Details	
Applications	Western Blot, Flow Cytometry, Flow (Intracellular), Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, CyTOF-ready
Recommended Dilutions	Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Flow (Intracellular), CyTOF-ready

Images

Flow Cytometry: CD3 zeta Antibody (H146-968) [PE] [NBP1-44683PE] - A cell surface stain was performed on hPBMCs with CD3 Zeta antibody (H146-968) NBP1-44683 (top image) and a matched isotype control (bottom image). Cells were incubated in an antibody dilution of 1:200 for 20 minutes at room temperature. Both antibodies were conjugated to PE.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-44683PE

NBP1-44683UV	CD3 zeta Antibody (H146-968) [DyLight 350]
NBP2-61462-20ug	Recombinant Human CD3 zeta His Protein
285-IF-100	IFN-gamma [Unconjugated]
NBL1-08925	CD3 zeta Overexpression Lysate

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-44683PE

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

