

Product Datasheet

LIME Antibody (LIME-06) - BSA Free NBP1-44945

Unit Size: 0.1 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-44945

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-44945



NBP1-44945

LIME Antibody (LIME-06) - BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	LIME-06
Preservative	15mM Sodium Azide
Isotype	IgG1 Kappa
Purity	Protein A purified
Buffer	Phosphate buffered saline (PBS), pH 7.4
Target Molecular Weight	30 kDa

Product Description	
Description	Novus Biologicals Mouse LIME Antibody (LIME-06) - BSA Free (NBP1-44945) is a monoclonal antibody validated for use in Flow and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	54923
Gene Symbol	LIME1
Species	Human
Specificity/Sensitivity	The antibody LIME-06 was raised against intracellular fragment corresponding to aa 141-295 of human LIME, a 30 kDa Lck-interacting transmembrane adaptor expressed by T cells.
Immunogen	Bacterially expressed intracellular fragment corresponding to aa 141-295 of human LIME.

Product Application Details	
Applications	Flow (Intracellular), Immunoprecipitation
Recommended Dilutions	Immunoprecipitation 1:50, Flow (Intracellular) 8-12 ug/ml

Images

Immunoprecipitation: LIME Antibody (LIME-06) [NBP1-44945] - Induction of LIME tyrosine phosphorylation. Peripheral blood T cells were left unstimulated (-) or stimulated with anti-human CD3 (MEM-92; cat. no. NB500-521) or anti-human CD4 (MEM-16), and LIME was immunoprecipitated from laurylmaltoside lysates with the LIME-06 antibody (immunoaffinity sorbent) and analyzed by Western blotting to visualize tyrosine-phosphorylated LIME (top) and total LIME (bottom).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-44945

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
202-IL-010	IL-2 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-44945

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

