

Product Datasheet

SRD5A2 Antibody NBP1-46510

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 3

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-46510

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-46510



NBP1-46510

SRD5A2 Antibody

Product Information	
Unit Size	0.1 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Reconstitution Instructions	Reconstitute in 0.1 ml of sterile water. Centrifuge to remove any insoluble material. Glycerol may be added (1:1) for additional stability. Please note the sample size is provided in reconstituted format.
Isotype	IgG
Purity	Unpurified
Buffer	Lyophilized from whole antisera


Product Description	
Description	Novus Biologicals Rabbit SRD5A2 Antibody (NBP1-46510) is a polyclonal antibody validated for use in IHC and WB. Anti-SRD5A2 Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	6716
Gene Symbol	SRD5A2
Species	Human, Mouse
Reactivity Notes	Other species not yet tested but expected to work in rat.
Immunogen	A synthetic peptide from amino acid region 190-250 of human SRD5A2 conjugated to blue carrier protein was used as the antigen.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1:4000, Immunohistochemistry 1:2000, Immunohistochemistry-Paraffin 1:2000
Application Notes	The optimal dilution should be determined by the end user. Although not tested this antibody may be useful in Immunohistochemistry-Paraffin/Frozen.




Images


Western Blot: SRD5A2 Antibody [NBP1-46510] - WB on mouse tissue lysate. Blocking: 1% LFDM for 30 min at RT; primary antibody: dilution 1:4000 incubated at 4C overnight.




Immunohistochemistry-Paraffin: SRD5A2 Antibody [NBP1-46510] - IHC-P on paraffin sections of human placenta. HIER: Tris-EDTA, pH 9 for 20 min. Blocking: 0.2% LFDM in TBST filtered thru 0.2 um. Detection was done using HRP polymer from manufacturers instructions; DAB chromogen. Primary antibody: dilution 1: 1000, incubated 30 min at RT using Autostainer. Sections were counterstained with Harris Hematoxylin.




Immunohistochemistry-Paraffin: SRD5A2 Antibody [NBP1-46510] - IHC-P on paraffin sections of human pancreas. HIER: Tris-EDTA, pH 9 for 20 min. Blocking: 0.2% LFDM in TBST filtered thru 0.2 um. Detection was done using HRP polymer following manufacturers instructions. DAB chromogen. Primary antibody: dilution 1: 1000, incubated 30 min at RT using Autostainer. Sections were counterstained with Harris Hematoxylin.




Immunohistochemistry-Paraffin: SRD5A2 Antibody [NBP1-46510] - IHC-P on paraffin sections of human pancreas. HIER: Tris-EDTA, pH 9 for 20 min. Blocking: 0.2% LFDM in TBST filtered thru 0.2 um. Detection was done using Novolink HRP polymer from Leica following manufacturers instructions; DAB chromogen. Primary antibody: dilution 1: 1000, incubated 30 min at RT using Autostainer. Sections were counterstained with Harris Hematoxylin.




Immunohistochemistry-Paraffin: SRD5A2 Antibody [NBP1-46510] - IHC-P on paraffin sections of human placenta. HIER: Tris-EDTA, pH 9 for 20 min. Blocking: 0.2% LFDM in TBST filtered thru 0.2 um. Detection was done using HRP polymer following manufacturers instructions; DAB chromogen. Primary antibody: dilution 1: 1000, incubated 30 min at RT using Autostainer. Sections were counterstained with Harris Hematoxylin.




Immunohistochemistry-Paraffin: SRD5A2 Antibody [NBP1-46510] - IHC-P on paraffin sections of human placenta. HIER: Tris-EDTA, pH 9 for 20 min. Blocking: 0.2% LFDM in TBST filtered thru 0.2 um. Detection was done using HRP polymer following manufacturers instructions; DAB chromogen. Primary antibody: dilution 1: 1000, incubated 30 min at RT using Autostainer. Sections were counterstained with Harris Hematoxylin.




Immunohistochemistry-Paraffin: SRD5A2 Antibody [NBP1-46510] - IHC-P on paraffin sections of human placenta. HIER: Tris-EDTA, pH 9 for 20 min. Blocking: 0.2% LFDM in TBST filtered thru 0.2 um. Detection was done using HRP polymer following manufacturers instructions; DAB chromogen. Primary antibody: dilution 1: 1000, incubated 30 min at RT using Autostainer. Sections were counterstained with Harris Hematoxylin.




Immunohistochemistry-Paraffin: SRD5A2 Antibody [NBP1-46510] - IHC-P on paraffin sections of human placenta. HIER: Tris-EDTA, pH 9 for 20 min. Blocking: 0.2% LFDM in TBST filtered thru 0.2 um. Detection was done using HRP polymer from manufacturers instructions; DAB chromogen. Primary antibody: dilution 1: 1000, incubated 30 min at RT using Autostainer. Sections were counterstained with Harris Hematoxylin.



Immunohistochemistry-Paraffin: SRD5A2 Antibody [NBP1-46510] - IHC-P on paraffin sections of human placenta. HIER: Tris-EDTA, pH 9 for 20 min. Blocking: 0.2% LFDM in TBST filtered thru 0.2 um. Detection was done using HRP polymer following manufacturers instructions; DAB chromogen. Primary antibody: dilution 1: 1000, incubated 30 min at RT using Autostainer. Sections were counterstained with Harris Hematoxylin.



Immunohistochemistry-Paraffin: SRD5A2 Antibody [NBP1-46510] - IHC-P on paraffin sections of human placenta. HIER: Tris-EDTA, pH 9 for 20 min. Blocking: 0.2% LFDM in TBST filtered thru 0.2 um. Detection was done using HRP polymer from manufacturers instructions; DAB chromogen. Primary antibody: dilution 1: 1000, incubated 30 min at RT using Autostainer. Sections were counterstained with Harris Hematoxylin.



Publications

Liu Z, Lin Z, Cao F et al. Upregulation of mir-1199-5p is associated with reduced type 2 5- α reductase expression in benign prostatic hyperplasia BMC urology 2022-11-07 [PMID: 36344974] (IHC-P, Human)

Details:

Dilution used in IHC 1:1500

Liu Z, Lin Z, Cao F et al. Upregulation of miR-1199-5p is associated with reduced type 2 5-alpha reductase in benign prostatic hyperplasia Research Square 2022-03-23 [PMID: 36344974] (IF/IHC, Human)

Lin ZM, Fan DD, Jin S et al. Methylated CpG dinucleotides in the 5-alpha reductase 2 gene may explain finasteride resistance in benign prostatic enlargement patients EBioMedicine 2020-11-23 [PMID: 33243958]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-46510

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-46510

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

