

Product Datasheet

MEK1 Antibody (OTI1F5) NBP1-47833

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-47833

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-47833



NBP1-47833


MEK1 Antibody (OT11F5)

Product Information	
Unit Size	0.1 ml
Concentration	0.99 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OT11F5
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.3), 1.0% BSA and 50% Glycerol
Target Molecular Weight	43.4 kDa
Product Description	
Description	Novus Biologicals Mouse MEK1 Antibody (OT11F5) (NBP1-47833) is a monoclonal antibody validated for use in IHC, WB and ICC/IF. Anti-MEK1 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	5604
Gene Symbol	MAP2K1
Species	Human, Mouse, Rat
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions.
Specificity/Sensitivity	This antibody is specific for Homo sapiens mitogen-activated protein kinase kinase 1 (MAP2K1).
Immunogen	Full-length protein expressed in 293T cell transfected with human MAP2K1 expression vector
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:2000, Immunohistochemistry 1:10-1:500, Immunocytochemistry/Immunofluorescence 1:100, Immunohistochemistry-Paraffin 1:50




Images


Western Blot: MEK1 Antibody (1F5) [NBP1-47833] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MEK1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MEK1.




Immunocytochemistry/Immunofluorescence: MEK1 Antibody (1F5) [NBP1-47833] Staining of COS7 cells transiently transfected by pCMV6-ENTRY MEK1.




Immunohistochemistry-Paraffin: MEK1 Antibody (1F5) [NBP1-47833] - Staining of paraffin-embedded prostate tissue using anti-MEK1 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: MEK1 Antibody (1F5) [NBP1-47833] - Staining of paraffin-embedded Adenocarcinoma of colon tissue using anti-MEK1 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: MEK1 Antibody (1F5) [NBP1-47833] - Staining of paraffin-embedded Carcinoma of bladder tissue using anti-MEK1 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: MEK1 Antibody (1F5) [NBP1-47833] - Staining of paraffin-embedded colon tissue using anti-MEK1 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: MEK1 Antibody (1F5) [NBP1-47833] - Staining of paraffin-embedded endometrium tissue using anti-MEK1 mouse monoclonal antibody.



Immunohistochemistry-Paraffin: MEK1 Antibody (1F5) [NBP1-47833] - Staining of paraffin-embedded Ovary tissue using anti-MEK1 mouse monoclonal antibody.



Immunohistochemistry-Paraffin: MEK1 Antibody (1F5) [NBP1-47833] - Staining of paraffin-embedded pancreas tissue using anti-MEK1 mouse monoclonal antibody.



Publications

Lu H Samanta D Xiang L et al. Chemotherapy triggers HIF-1-dependent glutathione synthesis and copper chelation that induces the breast cancer stem cell phenotype. Proc Natl Acad Sci U S A. 2015-07-30 [PMID: 26229077] (WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-47833

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-47833

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

