

Product Datasheet

CD8 Antibody (53-6.7) - BSA Free NBP1-49045

Unit Size: 0.5 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 3 **Publications: 28**

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-49045

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-49045



NBP1-49045

CD8 Antibody (53-6.7) - BSA Free

Product Information	
Unit Size	0.5 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	53-6.7
Preservative	0.02% Sodium Azide
Isotype	IgG2a Kappa
Purity	Protein A or G purified
Buffer	PBS
Target Molecular Weight	27 kDa

Product Description	
Description	This CD8 alpha antibody serves as an effective marker of cytotoxic T lymphocytes by binding to the CD8 co-receptor expressed on the cell surface of cytotoxic T cells, recognizing the topological domain of CD8 alpha. Because this CD8 antibody is made to the alpha chain it will recognize both the CD8 alpha - CD8 beta heterodimer, the most common form, as well as the CD8 alpha - CD8 alpha homodimer. The CD8A gene is also expressed in natural killer cells (NK cells), dendritic cells and cortical thymocytes making the CD8 alpha antibody a potential marker for these cells.
Host	Rat
Gene ID	925
Gene Symbol	CD8A
Species	Mouse, Rat
Immunogen	CD8 Antibody (53-6.7) was developed against mouse thymus or spleen.

Product Application Details	
Applications	Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation, Cell depletion, CyTOF-ready, Inhibition of T Cell Function
Recommended Dilutions	Flow Cytometry 1:10 - 1:1000, Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence 1:10-1:500. Use reported in scientific literature, Immunoprecipitation 1:10 - 1:500. Use reported in scientific literature (PMID 24565643), Immunohistochemistry-Paraffin 1:10-1:500, Immunohistochemistry-Frozen 1:10-1:500, Cell depletion, Inhibition of T Cell Function, CyTOF-ready




Application Notes


Each lot of this CD8a antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For immunofluorescent staining, the suggested use of this reagent is $<0.25 \text{ ug}/10^6$ cells in 100 μL volume. It is recommended that the reagent be titrated for optimal performance for each application. The 53-6.7 antibody has been reported to block antigen presentation via MHC class I and inhibit T cell responses to IL-2. This antibody has also been used for depletion of CD8a+ cells. Additional reported applications (for the relevant formats) include: immunoprecipitation, in vivo and in vitro cell depletion, inhibition of CD8 T cell proliferation, blocking of cytotoxicity, and immunohistochemical staining of both acetone-fixed frozen sections and zinc-fixed paraffin-embedded sections.

Images

Flow Cytometry: CD8 Antibody (53-6.7) [NBP1-49045] - CD8 alpha Antibody (53-6.7) [NBP1-49045] - Analysis of lymph nodes by multiple staining.




Immunohistochemistry-Paraffin: CD8 Antibody (53-6.7) [NBP1-49045] - CD8 alpha Antibody (53-6.7) [NBP1-49045] - CD8 alpha expression in mouse spleen tissue using anti-CD8 alpha antibody. Image from verified customer review.




Immunocytochemistry/Immunofluorescence: CD8 Antibody (53-6.7) [NBP1-49045] - CD8 alpha Antibody (53-6.7) [NBP1-49045] - Analysis of bone marrow tissue by single and multiple staining.




Flow Cytometry: CD8 Antibody (53-6.7) [NBP1-49045] - CD8 alpha Antibody (53-6.7) [NBP1-49045] - Staining of peripheral blood leukocytes.




Immunocytochemistry/Immunofluorescence: CD8 Antibody (53-6.7) [NBP1-49045] - CD8 alpha Antibody (53-6.7) [NBP1-49045] - Analysis of immersion fixed splenocytes. Primary antibody was used at a dilution of 10 ug/mL and incubated for 3 hours at room temperature.




Immunohistochemistry: CD8 Antibody (53-6.7) [NBP1-49045] - Histopathological analysis of spinal cord in EAE mice. Representative spinal cord sections from wild-type (WT), Asic1^{-/-}, and Asic2^{-/-} mice at 45 days after immunization. Sections of the spinal cord anterior horn were stained with hematoxylin and eosin (E&H) to assess inflammation and immunostained for MCHII, CD4⁺, CD8⁺ infiltrating cells. Image collected and cropped by CiteAb from the following publication (<https://onlinelibrary.wiley.com/doi/abs/10.1111/ejn.14302>), licensed under a CC-BY license.




Flow Cytometry: CD8 Antibody (53-6.7) - BSA Free [NBP1-49045] - Analysis of isolated mouse splenic CD8⁺ T cells using CD8 antibody (53-6.7) [Allophycocyanin] (NBP1-49045APC). Image from verified customer review.




Immunohistochemistry-Frozen: CD8 Antibody (53-6.7) [NBP1-49045] - CD8 alpha Antibody (53-6.7) [NBP1-49045] - Analysis of mouse spleen sections. CD8+ lymphocytes are marked by brown labeling of the cell surface.




Immunohistochemistry-Frozen: CD8 Antibody (53-6.7) [NBP1-49045] - Analysis of CD8+ T cells in allogeneic skin grafted onto a mouse.



Flow Cytometry: CD8 Antibody (53-6.7) [NBP1-49045] - Analysis of fixed murine splenocytes by multiple staining.



Immunohistochemistry-Paraffin: CD8 Antibody (53-6.7) [NBP1-49045] - CD8 alpha Antibody (53-6.7) [NBP1-49045] - CD8 alpha expression in human spleen tissue using anti-CD8 alpha antibody.



Publications

Chakrabarti J, Holokai L, Syu L et al. Hedgehog signaling induces PD-L1 expression and tumor cell proliferation in gastric cancer. *Oncotarget* 2018-12-21 [PMID: 30647844]

Won S, Lee C, Bae S et al. Mass-produced gram-negative bacterial outer membrane vesicles activate cancer antigen-specific stem-like CD8(+) T cells which enables an effective combination immunotherapy with anti-PD-1 *Journal of Extracellular Vesicles* 2023-08-10 [PMID: 37563797]

Zhou J, Wang H, Shu T et al. Myeloid-intrinsic cell cycle-related kinase drives immunosuppression to promote tumorigenesis *iScience* 2023-10-20 [PMID: 37731616]

Cai EP, Ishikawa Y, Zhang W et al. Genome-scale in vivo CRISPR screen identifies RNLS as a target for beta cell protection in type 1 diabetes *Nat Metab* 2020-07-27 [PMID: 32719542]

Mendes da Silva A Developing New Pharmacological Tools to Modulate Transposable Elements' Activity in Colorectal Cancer Thesis 2023-01-01 (ICC/IF)

Son S, Nam J, Kim AS et al. Induction of T-helper-17-cell-mediated anti-tumour immunity by pathogen-mimicking polymer nanoparticles *Nature biomedical engineering* 2022-12-23 [PMID: 36564626] (IHC-Fr, Mouse)

Salman S, Meyers DJ, Wicks EE Et al. HIF inhibitor 32-134D eradicates murine hepatocellular carcinoma in combination with anti-PD1 therapy *J Clin Invest* 2022-05-02 [PMID: 35499076] (FLOW, Mouse)

Details:

Citation using the Alexa Fluor 488 version of this antibody.

Bhatia S Protease activated nanosensors for the noninvasive diagnosis of community-acquired pneumonia Thesis 2022-01-01 (IHC-Fr, Mouse)

Reyes AF, Goldusky J, Bhimalli P et al. Tracking fluorescently labeled IL-15 and anti-PD-1 in the tumor microenvironment and draining lymph nodes *Journal of immunological methods* 2022-03-28 [PMID: 35358495] (IF/IHC, Mouse)

Jie X, Chen Y, Zhao Y et al. Targeting KDM4C enhances CD8+ T cell mediated antitumor immunity by activating chemokine CXCL10 transcription in lung cancer *Journal for immunotherapy of cancer* 2022-02-01 [PMID: 35121645] (IF/IHC, Mouse)

Zhang Z, Na H, Gan Q Et al. Monomeric C-reactive protein via endothelial CD31 for neurovascular inflammation in an ApoE genotype-dependent pattern: A risk factor for Alzheimer's disease? *Aging cell* 2021-10-23 [PMID: 34687487] (IHC-Fr, Mouse)

Yang W, Feng Y, Zhou J et al. A selective HDAC8 inhibitor potentiates antitumor immunity and efficacy of immune checkpoint blockade in hepatocellular carcinoma *Science translational medicine* 2021-04-07 [PMID: 33827976]

More publications at <http://www.novusbio.com/NBP1-49045>





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-49045

HAF005	Goat anti-Rat IgG Secondary Antibody [HRP]
NB7115	Goat anti-Rat IgG (H+L) Secondary Antibody [HRP]
NBP1-43321-0.5mg	Rat IgG2a Kappa Light Chain Isotype Control (R2a)
NBP1-49045PE	CD8 Antibody (53-6.7) [PE]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-49045

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

