

# Product Datasheet

## IGF2BP3 Antibody (8F11) NBP1-51674

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-51674](http://www.novusbio.com/NBP1-51674)

Updated 1/11/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-51674](http://www.novusbio.com/reviews/destination/NBP1-51674)



**NBP1-51674****IGF2BP3 Antibody (8F11)**

| <b>Product Information</b>     |                                                                                        |
|--------------------------------|----------------------------------------------------------------------------------------|
| <b>Unit Size</b>               | 0.1 ml                                                                                 |
| <b>Concentration</b>           | This product is unpurified. The exact concentration of antibody is not quantifiable.   |
| <b>Storage</b>                 | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| <b>Clonality</b>               | Monoclonal                                                                             |
| <b>Clone</b>                   | 8F11                                                                                   |
| <b>Preservative</b>            | 0.03% Sodium Azide                                                                     |
| <b>Isotype</b>                 | IgG1                                                                                   |
| <b>Purity</b>                  | Ascites                                                                                |
| <b>Buffer</b>                  | Ascites                                                                                |
| <b>Target Molecular Weight</b> | 69 kDa                                                                                 |


| <b>Product Description</b> |                                                                                                                                                                                                   |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Novus Biologicals Mouse IGF2BP3 Antibody (8F11) (NBP1-51674) is a monoclonal antibody validated for use in IHC, WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| <b>Host</b>                | Mouse                                                                                                                                                                                             |
| <b>Gene ID</b>             | 10643                                                                                                                                                                                             |
| <b>Gene Symbol</b>         | IGF2BP3                                                                                                                                                                                           |
| <b>Species</b>             | Human                                                                                                                                                                                             |
| <b>Immunogen</b>           | Purified recombinant fragment of human IGF2BP3 expressed in E. Coli.                                                                                                                              |

| <b>Product Application Details</b> |                                                                                                                         |
|------------------------------------|-------------------------------------------------------------------------------------------------------------------------|
| <b>Applications</b>                | Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunohistochemistry                                                |
| <b>Recommended Dilutions</b>       | Western Blot 1:500-1:2000, ELISA 1:10000, Immunohistochemistry 1:200-1:1000, Immunohistochemistry-Paraffin 1:200-1:1000 |




## Images

Western Blot: IGF2BP3 Antibody (8F11) [NBP1-51674] - Analysis using IGF2BP3 mouse mAb against JURKET (1), K562 (2) and NTERA-2 (3) cell lysate.



Immunohistochemistry-Paraffin: IGF2BP3 Antibody (8F11) [NBP1-51674] - Analysis of cancer (left) and colon tumour tissues (right) using IGC2BP3 mouse mAb with DAB staining.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-51674**

---

|                  |                                                    |
|------------------|----------------------------------------------------|
| NBP2-33376H      | Blue Marker Antibody (6F4-F6) [HRP]                |
| HAF007           | Goat anti-Mouse IgG Secondary Antibody [HRP]       |
| NB7539           | Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP] |
| NBP1-97005-0.5mg | Mouse IgG1 Isotype Control (MG1)                   |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-51674](http://www.novusbio.com/reviews/submit/NBP1-51674)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

