

# Product Datasheet

## RRP9 Antibody - BSA Free NBP1-57217

Unit Size: 100ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-57217](http://www.novusbio.com/NBP1-57217)

Updated 2/21/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-57217](http://www.novusbio.com/reviews/destination/NBP1-57217)



**NBP1-57217**

RRP9 Antibody - BSA Free

**Product Information**

|                      |                                                                                        |
|----------------------|----------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 100ul                                                                                  |
| <b>Concentration</b> | 1 mg/ml                                                                                |
| <b>Storage</b>       | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| <b>Clonality</b>     | Polyclonal                                                                             |
| <b>Preservative</b>  | 0.09% Sodium Azide                                                                     |
| <b>Isotype</b>       | IgG                                                                                    |
| <b>Purity</b>        | Protein A purified                                                                     |
| <b>Buffer</b>        | PBS, 2% Sucrose                                                                        |

**Product Description**


|                    |                                                                                                                                                                                                                                                                                                                                      |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | The addition of 50% glycerol is optional for those storing this antibody at -20C and not aliquoting smaller units. However, please note that glycerol may interrupt some downstream antibody applications and should be added with caution.                                                                                          |
| <b>Host</b>        | Rabbit                                                                                                                                                                                                                                                                                                                               |
| <b>Gene ID</b>     | 9136                                                                                                                                                                                                                                                                                                                                 |
| <b>Gene Symbol</b> | RRP9                                                                                                                                                                                                                                                                                                                                 |
| <b>Species</b>     | Human                                                                                                                                                                                                                                                                                                                                |
| <b>Immunogen</b>   | Synthetic peptides corresponding to RRP9 (RRP9, small subunit (SSU) processome component, homolog (yeast)) The peptide sequence was selected from the middle region of RRP9. Peptide sequence IPRAKKGAEGKPPGHSSHVLCMAISSDGKYLASGDRSKLILIWEAQSCQH The peptide sequence for this immunogen was taken from within the described region. |

**Product Application Details**

|                              |                        |
|------------------------------|------------------------|
| <b>Applications</b>          | Western Blot           |
| <b>Recommended Dilutions</b> | Western Blot 1.0 ug/ml |

**Images**

Western Blot: RRP9 Antibody [NBP1-57217] - Jurkat cell lysate, concentration 5.0ug/ml.



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP1-57217**

---

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891         | Rabbit IgG Isotype Control                          |
| H00009136-P01-10ug | Recombinant Human RRP9 GST (N-Term) Protein         |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-57217](http://www.novusbio.com/reviews/submit/NBP1-57217)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

