

Product Datasheet

DAPK3/ZIPK Antibody - BSA Free NBP1-76962

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-76962

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-76962



NBP1-76962


DAPK3/ZIPK Antibody - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	27 kDa
Product Description	
Description	Novus Biologicals Rabbit DAPK3/ZIPK Antibody - BSA Free (NBP1-76962) is a polyclonal antibody validated for use in IHC, WB, ELISA and ICC/IF. Anti-DAPK3/ZIPK Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	1613
Gene Symbol	DAPK3
Species	Human, Mouse, Rat
Immunogen	Antibody was raised against a peptide corresponding to amino acids near the center of human ZIP kinase. The immunogen is located within amino acids 270 - 320 of ZIPK. Amino Acid Sequence: RRRNVRGEDSGRKPERRRLK
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1 ug/ml, ELISA 1:100-1:2000, Immunocytochemistry/ Immunofluorescence 5 ug/mL
Application Notes	An approximate 40 kDa band can be detected in Western Blot, which represents the pro-enzyme of DNase II.




Images


Western Blot: DAPK3/ZIPK Antibody [NBP1-76962] - Analysis of ZIP kinase in (A) HeLa and (B) Jurkat lysates with ZIP kinase antibody at 1 ug/mL.



Immunocytochemistry/Immunofluorescence: DAPK3/ZIPK Antibody [NBP1-76962] - Staining of ZIPK in Jurkat cells with ZIPK antibody at 5 ug/ml.



Immunocytochemistry/ Immunofluorescence: DAPK3/ZIPK Antibody - BSA Free [NBP1-76962] - Immunofluorescence of DAPK3/ZIPK in Jurkat cells with DAPK3/ZIPK antibody at 10 ug/mL.



Publications

Chen HM, MacDonald JA Death-associated protein kinases and intestinal epithelial homeostasis *Anatomical record* (Hoboken, N.J. : 2007) 2022-06-23 [PMID: 35735750] (IHC-P)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-76962

NBP1-76962PEP	DAPK3/ZIPK Antibody Blocking Peptide
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-76962

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

