

# Product Datasheet

## TetR Antibody [FITC] NBP1-78424F

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-78424F](http://www.novusbio.com/NBP1-78424F)

Updated 10/29/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-78424F](http://www.novusbio.com/reviews/destination/NBP1-78424F)



**NBP1-78424F**

TetR Antibody [FITC]

| Product Information          |   |
|------------------------------|---|
| <b>Unit Size</b>             | 0.1 ml  |
| <b>Concentration</b>         | Please see the vial label for concentration. If unlisted please contact technical services.                                 |
| <b>Storage</b>               | Store at 4C in the dark.  |
| <b>Clonality</b>             | Polyclonal  |
| <b>Preservative</b>          | 0.05% Sodium Azide  |
| <b>Isotype</b>               | IgG   |
| <b>Conjugate</b>             | FITC  |
| <b>Purity</b>                | Immunogen affinity purified   |
| <b>Buffer</b>                | PBS   |
| Product Description          |   |
| <b>Host</b>                  | Rabbit  |
| <b>Gene ID</b>               | 4924774   |
| <b>Species</b>               | Non-species specific  |
| <b>Reactivity Notes</b>      | This antibody is used to monitor the regulation of gene expression in bacteria cells.                                       |
| <b>Immunogen</b>             | A synthetic peptide made to an N-terminal portion of Escherichia coli TetR protein (between residues 1-50) [UniProt P04483] |
| Product Application Details  |   |
| <b>Applications</b>          | Western Blot  |
| <b>Recommended Dilutions</b> | Western Blot  |
| <b>Application Notes</b>     | This TetR antibody is useful for Western blot.  |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-78424F**

---

|               |                                   |
|---------------|-----------------------------------|
| NBP2-24892    | Rabbit IgG Isotype Control [FITC] |
| NB600-234PEP  | TetR Antibody Blocking Peptide    |
| NBP1-78424PEP | TetR Antibody Blocking Peptide    |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-78424F](http://www.novusbio.com/reviews/submit/NBP1-78424F)

Earn gift cards/discounts by submitting a publication using this product:

[www.novusbio.com/publications](http://www.novusbio.com/publications)

