

Product Datasheet

SLFN5 Antibody - BSA Free

NBP1-81178

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-81178

Updated 3/4/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-81178



NBP1-81178

SLFN5 Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description

Host	Rabbit
Gene ID	162394
Gene Symbol	SLFN5
Species	Human
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: ERHGVGLDVPPIFRSHLDKMQKENHFLIFVKSWNTEAGVPLATLCSNLYHRER TSTDVMDSQEALAFKCRQTPTNINVSNSLGPQAAQGSVQYEGNINVSAAL FDRKRLQYLEKLNLPETHVEFVMFSTDVSHCVKD


Product Application Details

Applications	Western Blot, Simple Western, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04 - 0.4 ug/ml, Simple Western 1:20, Immunohistochemistry 1:500 - 1:1000, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:500 - 1:1000
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100. In Simple Western only 10 - 15 uL of the recommended dilution is used per data point. See Simple Western Antibody Database for Simple Western validation: Tested in RT-4 and U-251MG, separated by Size, antibody dilution of 1:20, apparent MW was 92 kDa. Separated by Size-Wes, Sally Sue/Peggy Sue.



Images


Western Blot: SLFN5 Antibody [NBP1-81178] - Analysis in human cell lines SK-MEL-30 and Caco-2. Corresponding RNA-seq data are presented for the same cell lines. Loading control: Anti-HDAC1.




Immunohistochemistry-Paraffin: SLFN5 Antibody [NBP1-81178] - Staining of human skeletal muscle shows weak nuclear positivity in myocytes as expected.




Immunohistochemistry-Paraffin: SLFN5 Antibody [NBP1-81178] - Staining of human Fallopian tube shows strong nuclear positivity in glandular cells.




Immunohistochemistry-Paraffin: SLFN5 Antibody [NBP1-81178] - Staining of human lymph node shows strong nuclear positivity in non-germinal center cells.




Immunohistochemistry-Paraffin: SLFN5 Antibody [NBP1-81178] - Staining of human endometrium shows strong nuclear positivity in glandular cells.A



Simple Western: SLFN5 Antibody [NBP1-81178] - Simple Western lane view shows a specific band for SLFN5 in 0.2 mg/ml of RT-4 (left) and U-251MG sp (right) lysate. This experiment was performed under reducing conditions using the 12-230 kDa separation system.



Simple Western: SLFN5 Antibody [NBP1-81178] - Electropherogram image(s) of corresponding Simple Western lane view. SLFN5 antibody was used at 1:20 dilution on RT-4 and U-251MG sp lysates(s).



Publications

Huang J, Wu C, Kloeber JA et al. SLFN5-mediated chromatin dynamics sculpt higher-order DNA repair topology Molecular cell 2023-02-25 [PMID: 36854302] (WB, ICC/IF, Human)

Katsoulidis E, Mavrommatis E, Woodard J et al. Role of interferon {alpha} (IFN{alpha})-inducible Schlafen-5 in regulation of anchorage-independent growth and invasion of malignant melanoma cells. J Biol Chem 2010-12-01 [PMID: 20956525]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-81178

NBP1-81178PEP	SLFN5 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-81178

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

