

Product Datasheet

NUAK2 Antibody - BSA Free

NBP1-81880

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 6

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-81880

Updated 12/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-81880



NBP1-81880

NUAK2 Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol


Product Description	
Description	Novus Biologicals Rabbit NUAK2 Antibody - BSA Free (NBP1-81880) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-NUAK2 Antibody: Cited in 5 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	81788
Gene Symbol	NUAK2
Species	Human, Mouse
Reactivity Notes	Rat (81%). Use in Mouse reported in scientific literature (PMID:31350328).
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: DPKEQKPPQASGLLLHRKGILKLNKFSQTALELAAPTTFGSLDELAPPRPLAR ASRPSGAVSEDSILSSESFDQLDLPERLPEPPLRGCVSVDNLTGLEEPPSEGP GSCLRRWRQDPLGDSCFSLTDCQEVATATYRQALRVCSKLT

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Knockdown Validated
Recommended Dilutions	Western Blot 0.04 - 0.4 ug/ml, Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:50-1:200, Knockdown Validated
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100.




Images


Western Blot: NUAK2 Antibody [NBP1-81880] - NUAK2 activity is involved in serum-induced or LPA-induced YAP/TAZ dephosphorylation and activation. FBS and LPA activate YAP/TAZ and induce NUAK2 expression. YAP/TAZ phosphorylation status was monitored using Phos-Tag gels, while NUAK2 and ANKRD1 mRNA expression was determined by qPCR at the indicated times. Expression data are plotted as the mean \pm SD ($n = 3$). Image collected and cropped by CiteAb from the following publication (<https://www.nature.com/articles/s41467-018-05939-2>), licensed under a CC-BY license.




Immunocytochemistry/Immunofluorescence: NUAK2 Antibody [NBP1-81880] - Staining of human cell line U-251 MG shows localization to nucleoplasm, nucleoli fibrillar center & cytosol. Antibody staining is shown in green.



Immunohistochemistry-Paraffin: NUAK2 Antibody [NBP1-81880] - Staining of human testis shows moderate to strong cytoplasmic and nuclear positivity in cells in seminiferous ducts.




Western Blot: NUAK2 Antibody [NBP1-81880] - Analysis in control (vector only transfected HEK293T lysate) and NUAK2 over-expression lysate (Co-expressed with a C-terminal myc-DDK tag (3.1 kDa) in mammalian HEK293T cells).




Immunohistochemistry-Paraffin: NUAK2 Antibody [NBP1-81880] - Staining of human cerebellum shows strong membranous positivity in Purkinje cells.




Immunohistochemistry-Paraffin: NUAK2 Antibody [NBP1-81880] - Staining of human cervix shows strong nuclear positivity in squamous epithelial cells.




Immunohistochemistry-Paraffin: NUAK2 Antibody [NBP1-81880] - Staining of human duodenum shows moderate to strong nuclear positivity in glandular cells.




Western Blot: NUAK2 Antibody [NBP1-81880] - NUAK2 regulates YAP/TAZ activity through LATS. NUAK2 requires LATS to regulate YAP/TAZ localization, phosphorylation, and target gene expression in MDA-MB231 cells. YAP/TAZ phosphorylation was monitored using Phos-Tag gels. Relative phosphorylation levels from blots is quantitated (right). Image collected and cropped by CiteAb from the following publication (<https://www.nature.com/articles/s41467-018-05939-2>), licensed under a CC-BY license.




Western Blot: NUAK2 Antibody [NBP1-81880] - NUAK2 activity is involved in serum-induced or LPA-induced YAP/TAZ dephosphorylation & activation. a, b FBS & LPA activate YAP/TAZ & induce NUAK2 expression. YAP/TAZ phosphorylation status was monitored using Phos-Tag gels (a), while NUAK2 & ANKRD1 mRNA expression was determined by qPCR at the indicated times. Expression data are plotted as the mean \pm SD ($n = 3$). c–h WZ4003 (10 μ M) blocks serum/LPA-induced NUAK2 & ANKRD1 mRNA expression. FBS-induced or LPA-induced YAP/TAZ dephosphorylation (e & g), or NUAK2 & ANKRD1 mRNA expression (c, f, h) was monitored at the indicated times (d). Expression data are plotted as the mean \pm range for a representative experiment (c, f, h) Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/30158528>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.




Western Blot: NUAK2 Antibody [NBP1-81880] - NUAK2 regulates YAP/TAZ activity through LATS. a Loss of NUAK2 enhances YAP/TAZ phosphorylation in MDA-MB231 cells. Quantitation of relative phosphorylation levels from blots is shown. b NUAK2 interacts with LATS1 & LATS2, but not MST1, MST2, or TAZ in the LUMIER protein interaction screen. c Lysates from HEK293T cells, transfected with NUAK2-Flag, were subjected to immunoprecipitation (IP) using anti-LATS1 antibody for endogenous LATS1 & co-immunoprecipitated NUAK2-Flag was detected by immunoblotting. Protein expression levels were confirmed (Totals). WT wild type, KR K81R, kinase dead. d–f NUAK2 requires LATS to regulate YAP/TAZ localization, phosphorylation, & target gene expression in MDA-MB231 cells. d YAP/TAZ localization was quantitated & plotted as the mean \pm SD ($n = 3$) with representative images shown on the right. Scale bars, 25 μ m. N: nuclear, C: cytoplasmic. e Target gene expression was determined by qPCR. Data are plotted as the mean \pm SD ($n = 3$). f YAP/TAZ phosphorylation was monitored using Phos-Tag gels. Relative phosphorylation levels from blots is quantitated (right) Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/30158528>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.




Western Blot: NUAK2 Antibody [NBP1-81880] - NUAK2 regulates YAP/TAZ activity through LATS. a Loss of NUAK2 enhances YAP/TAZ phosphorylation in MDA-MB231 cells. Quantitation of relative phosphorylation levels from blots is shown. b NUAK2 interacts with LATS1 & LATS2, but not MST1, MST2, or TAZ in the LUMIER protein interaction screen. c Lysates from HEK293T cells, transfected with NUAK2-Flag, were subjected to immunoprecipitation (IP) using anti-LATS1 antibody for endogenous LATS1 & co-immunoprecipitated NUAK2-Flag was detected by immunoblotting. Protein expression levels were confirmed (Totals). WT wild type, KR K81R, kinase dead. d–f NUAK2 requires LATS to regulate YAP/TAZ localization, phosphorylation, & target gene expression in MDA-MB231 cells. d YAP/TAZ localization was quantitated & plotted as the mean \pm SD ($n = 3$) with representative images shown on the right. Scale bars, 25 μ m. N: nuclear, C: cytoplasmic. e Target gene expression was determined by qPCR. Data are plotted as the mean \pm SD ($n = 3$). f YAP/TAZ phosphorylation was monitored using Phos-Tag gels. Relative phosphorylation levels from blots is quantitated (right) Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/30158528>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Western Blot: NUAK2 Antibody [NBP1-81880] - NUAK2 is associated with high-grade human bladder cancer. a, c Expression of NUAK2 & the YAP Signature Genes (YAP SG) is enhanced both in high-grade non-muscle-invasive (HG-NMIBC; n = 13, **p = 0.0071) & muscle-invasive (pT2 MIBC; n = 9, **p = 0.0033) bladder cancers³⁹ as compared to LG-NMIBC (n = 27) using an unpaired t test. b Elevated expression of NUAK2 is correlated with disease recurrence in a cohort⁴¹ of MIBC patient samples with n = 57 disease free & n = 56, recurred samples (*p = 0.0172). Box & whisker plots show the median (line in the box), first & third quartiles (lower & upper ends of the box), & the minimum & maximum values (whiskers in the plot). Dots represent a single patient sample. d A heat map depicting expression of NUAK2 & YAP SG in NMIBC LG (open circle), HG (closed circle), & MIBC pT2 (slashed circle) bladder cancer samples³⁹. e Characterization of high-grade-derived & low-grade-derived bladder cancer cell lines. Scale bar, 25 μ m. f NUAK2 mRNA expression is elevated in HG-derived (TCCSUP & T24) BC cell lines as compared to LG lines. Data are plotted as the mean \pm SD (n = 3). g–i Blocking NUAK2 activity with 10 μ M WZ4003 (g) or NUAK2 expression in stable clones (#5 or #7) using inducible NUAK2 shRNAs (h, i) inhibits the growth of HG BC cell lines as measured by SRB (g) or DAPI (h) staining. Data are plotted as the mean \pm SD of a representative experiment, n = 3 (for TCC) or mean \pm range of two experiments (for T24) Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/30158528>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Western Blot: NUAK2 Antibody [NBP1-81880] - NUAK2 regulates YAP/TAZ activity through LATS. a Loss of NUAK2 enhances YAP/TAZ phosphorylation in MDA-MB231 cells. Quantitation of relative phosphorylation levels from blots is shown. b NUAK2 interacts with LATS1 & LATS2, but not MST1, MST2, or TAZ in the LUMIER protein interaction screen. c Lysates from HEK293T cells, transfected with NUAK2-Flag, were subjected to immunoprecipitation (IP) using anti-LATS1 antibody for endogenous LATS1 & co-immunoprecipitated NUAK2-Flag was detected by immunoblotting. Protein expression levels were confirmed (Totals). WT wild type, KR K81R, kinase dead. d–f NUAK2 requires LATS to regulate YAP/TAZ localization, phosphorylation, & target gene expression in MDA-MB231 cells. d YAP/TAZ localization was quantitated & plotted as the mean \pm SD (n = 3) with representative images shown on the right. Scale bars, 25 μ m. N: nuclear, C: cytoplasmic. e Target gene expression was determined by qPCR. Data are plotted as the mean \pm SD (n = 3). f YAP/TAZ phosphorylation was monitored using Phos-Tag gels. Relative phosphorylation levels from blots is quantitated (right) Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/30158528>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Publications

Jo H, Munoz S, Dalvi A et al. A fetal oncogene NUA2 is an emerging therapeutic target in glioblastoma EMBO Molecular Medicine 2025-08-06 [PMID: 40770117]

Fu W, Zhao MT, Driver LM Et al. NUA family kinase 2 is a novel therapeutic target for prostate cancer Molecular carcinogenesis 2021-11-24 [PMID: 34818445] (WB, Human)

Hollstein P E, Eichner L J et al. The AMPK-Related Kinases SIK1 and SIK3 Mediate Key Tumor-Suppressive Effects of LKB1 in NSCLC. Cancer Discov 2019-01-11 [PMID: 31350328] (WB, Mouse)

Gill MK, Christova T, Zhang YY et al. A feed forward loop enforces YAP/TAZ signaling during tumorigenesis Nat Commun 2018-08-29 [PMID: 30158528] (WB, Human)

Stadler C, Rexhepaj E, Singan VR et al. Immunofluorescence and fluorescent-protein tagging show high correlation for protein localization in mammalian cells. Nat Methods 2013-04-01 [PMID: 23435261]

Kim JG, Lee SJ, Chae YS et al. Association between Phosphorylated AMP-Activated Protein Kinase and MAPK3/1 Expression and Prognosis for Patients with Gastric Cancer. Oncology 2013-07-16 [PMID: 23860205]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-81880

NBP1-81880PEP	NUAK2 Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-81880

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

