

Product Datasheet

PSAP Antibody - BSA Free NBP1-83363

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-83363

Updated 12/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-83363



NBP1-83363

PSAP Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Description	Novus Biologicals Rabbit PSAP Antibody - BSA Free (NBP1-83363) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5660
Gene Symbol	PSAP
Species	Human
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: TNSTFVQALVEHVKEECDRLGPGMADICKNYISQYSEIAIQMMMHHMQPKEICAL VGFCDEVKEMPMQTLVPAKVASKNVIPALELVEPIKKHEVPAKSDVYCEVCEFL VKEVTKLIDNNKTEKEILDAFDKMCSKLPKSLSEE

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:200 - 1:500, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:200 - 1:500
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF, Fixation Permeabilization: Use PFA/Triton X-100.

Images


Immunohistochemistry-Paraffin: PSAP Antibody [NBP1-83363] - Staining in human cerebral cortex and skeletal muscle tissues using NBP1-83363 antibody. Corresponding PSAP RNA-seq data are presented for the same tissues.




Immunocytochemistry/Immunofluorescence: PSAP Antibody [NBP1-83363] - Staining of human cell line U-2 OS shows localization to vesicles. Antibody staining is shown in green.



Immunohistochemistry-Paraffin: PSAP Antibody [NBP1-83363] - Staining of human cerebral cortex shows strong granular cytoplasmic positivity in neurons.




Immunohistochemistry-Paraffin: PSAP Antibody [NBP1-83363] - Staining of human placenta shows strong granular cytoplasmic positivity in trophoblastic cells.




Immunohistochemistry-Paraffin: PSAP Antibody [NBP1-83363] - Staining of human prostate shows strong granular cytoplasmic positivity in glandular cells.



Immunohistochemistry-Paraffin: PSAP Antibody [NBP1-83363] - Staining of human skeletal muscle shows very weak positivity in myocytes as expected.



Analysis in human cell lines SK-MEL-30 and HEK293 using Anti-PSAP antibody. Corresponding PSAP RNA-seq data are presented for the same cell lines. Loading control: Anti-PFN1.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-83363

NBP1-83363PEP	PSAP Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-83363

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

