

Product Datasheet

SNRK Antibody - BSA Free

NBP1-83668

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-83668

Updated 12/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-83668



NBP1-83668

SNRK Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Description	Novus Biologicals Rabbit SNRK Antibody - BSA Free (NBP1-83668) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	54861
Gene Symbol	SNRK
Species	Human, Mouse, Rat
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: VPQDLEDDLTTATPLSHATVPQSPARAADSVLNGHRSKGLCDSAKKDDLPELAG PALSTVPPASLKPTASGRKCLFRVEEDEDEEEDKKPMSLSTQVVLRRKPSVT NRLTSRKSAVLNQIFEEGESDDEFDMDENLPPKLSRLKM

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:50 - 1:200
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization, Use PFA/Triton X-100.

Images


Immunocytochemistry/Immunofluorescence: SNRK Antibody [NBP1-83668] - Analysis of human brain microvascular endothelial cells using SNRK antibody. SNRK overexpression is found to be highly expressed in the nucleus, SNRK RED, DAPI BLUE. Image from verified customer review.




Immunocytochemistry/Immunofluorescence: SNRK Antibody [NBP1-83668] - Staining of human cell line U-2 OS shows localization to cell junctions & vesicles. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: SNRK Antibody [NBP1-83668] -Staining of human pancreas shows strong nuclear positivity in exocrine glandular cells and endocrine glandular cells.




Immunohistochemistry-Paraffin: SNRK Antibody [NBP1-83668] -Staining of human lymphoid tissues shows strong nuclear positivity in non-germinal center cells.




Immunohistochemistry-Paraffin: SNRK Antibody [NBP1-83668] -Staining of human testis shows strong nucleus positivity in endothelial cells, leydig cells and moderate nucleus positivity in cells in seminiferous cells.




Immunohistochemistry-Paraffin: SNRK Antibody [NBP1-83668] -Staining of human skeletal muscle shows strong nuclear positivity in myocytes.



Analysis in human cell lines PC-3 and MCF-7 using Anti-SNRK antibody. Corresponding SNRK RNA-seq data are presented for the same cell lines. Loading control: Anti-PPIB.



Analysis in mouse cell line NIH-3T3 and rat cell line NBT-II.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-83668

NBP1-83668PEP	SNRK Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-83668

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

