

Product Datasheet

ATAD2 Antibody - BSA Free

NBP1-84123

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-84123

Updated 12/2/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-84123



NBP1-84123

ATAD2 Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol


Product Description	
Description	Novus Biologicals Rabbit ATAD2 Antibody - BSA Free (NBP1-84123) is a polyclonal antibody validated for use in IHC, WB, ICC/IF and Simple Western. Anti-ATAD2 Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	29028
Gene Symbol	ATAD2
Species	Human
Reactivity Notes	Immunogen displays the following percentage of sequence identity for non-tested species: Mouse (82%), Rat (82%). Reactivity reported in scientific literature (PMID: 23393560)
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: TAYAIKEELDEDFEQLCEEIQESRKKRGCSSSKYAPSYHVMPKQNSTLVGDK RSDPEQNEKLKTPSTPVACSTPAQLKRKIRKKS NWYLGTIKKRRKISQAKDDSQ NAIDHK

Product Application Details	
Applications	Western Blot, Simple Western, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Knockdown Validated
Recommended Dilutions	Western Blot 0.4 ug/ml, Simple Western 1:50, Immunohistochemistry 1:500 - 1:1000, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:500-1:1000, Knockdown Validated
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100. In Simple Western only 10 - 15 uL of the recommended dilution is used per data point. See Simple Western Antibody Database for Simple Western validation: Tested in U-251MG sp, separated by Size, antibody dilution of 1:50, apparent MW was 178 kDa. Separated by Size-Wes, Sally Sue/Peggy Sue.




Images


Immunohistochemistry-Paraffin: ATAD2 Antibody [NBP1-84123] - Analysis in human testis and skeletal muscle tissues using NBP1-84123 antibody. Corresponding ATAD2 RNA-seq data are presented for the same tissues.




Immunocytochemistry/Immunofluorescence: ATAD2 Antibody [NBP1-84123] - Staining of human cell line U-2 OS shows localization to nucleoplasm. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: ATAD2 Antibody [NBP1-84123] - Staining of human Skeletal muscle shows very weak positivity in myocytes.




Immunohistochemistry-Paraffin: ATAD2 Antibody [NBP1-84123] - Staining of human Bone marrow shows strong nuclear positivity in hematopoietic cells.




Immunohistochemistry-Paraffin: ATAD2 Antibody [NBP1-84123] - Staining of human Lymph node shows moderate nuclear positivity in non-germinal center cells.




Immunohistochemistry-Paraffin: ATAD2 Antibody [NBP1-84123] - Staining of human Testis shows strong nuclear positivity in cells in seminiferous ducts.



Simple Western: ATAD2 Antibody [NBP1-84123] - Simple Western lane view shows a specific band for ATAD2 in 0.2 mg/ml of U-251MG sp lysate. This experiment was performed under reducing conditions using the 66-440 kDa separation system.



Analysis in U2OS cells) transfected with control siRNA, target specific siRNA probe #1 and #2, using Anti-ATAD2 antibody. Remaining relative intensity is presented.



Publications

Wan WN, Zhang YX, Wang XM et al. ATAD2 is highly expressed in ovarian carcinomas and indicates poor prognosis. *Asian Pac J Cancer Prev* 2014-01-01 [PMID: 24761900] (IF/IHC, Human)

Raeder MB, Birkeland E, Trovik J et al. Integrated Genomic Analysis of the 8q24 Amplification in Endometrial Cancers Identifies ATAD2 as Essential to MYC-Dependent Cancers. *PLoS One* 2013-01-01 [PMID: 23393560]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to NBP1-84123

NBP1-84123PEP	ATAD2 Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-84123

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



