

Product Datasheet

ARHGAP30 Antibody - BSA Free

NBP1-84653

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-84653

Updated 12/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-84653



NBP1-84653


ARHGAP30 Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol
Product Description	
Description	Novus Biologicals Rabbit ARHGAP30 Antibody - BSA Free (NBP1-84653) is a polyclonal antibody validated for use in IHC and ICC/IF. Anti-ARHGAP30 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	257106
Gene Symbol	ARHGAP30
Species	Human
Reactivity Notes	Immunogen displays the following percentage of sequence identity for non-tested species: Mouse (81%), Rat (81%)
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: ANPCPRPGRLDGTGPGERAWGSRASRSSWRNGGSLSFDAVALARDRQRTEA QGVRRRTQTCTEGGDYCLIPRTSPCSMISAHSPRPLSCLELPSEGAEGSGSRSR
Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Immunohistochemistry 1:200 - 1:500, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:200 - 1:500
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100.




Images


Immunocytochemistry/Immunofluorescence: ARHGAP30 Antibody [NBP1-84653] - Staining of human cell line A-431 shows localization to vesicles. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: ARHGAP30 Antibody [NBP1-84653] - Staining of human skeletal muscle shows no positivity in myocytes as expected.




Immunohistochemistry-Paraffin: ARHGAP30 Antibody [NBP1-84653] - Staining of human tonsil shows strong cytoplasmic positivity in non-germinal center cells.




Immunohistochemistry-Paraffin: ARHGAP30 Antibody [NBP1-84653] - Staining of human lung shows moderate cytoplasmic positivity in macrophages.



Immunohistochemistry-Paraffin: ARHGAP30 Antibody [NBP1-84653] - Staining of human pancreas shows no positivity in exocrine glandular cells as expected.



Analysis in human tonsil and skeletal muscle tissues using HPA036300 antibody. Corresponding ARHGAP30 RNA-seq data are presented for the same tissues.



Publications

Zhou Y, Hua Z, Zhu Y et al. Upregulation of ARHGAP30 attenuates pancreatic cancer progression by inactivating the beta-catenin pathway Cancer Cell Int. 2020-06-10 [PMID: 32536813] (IHC-P, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-84653

NBP1-84653PEP	ARHGAP30 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-84653

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

