

Product Datasheet

Dynactin Subunit 2/DCTN2/DCTN-50 Antibody - BSA Free NBP1-85277

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-85277

Updated 3/4/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-85277



NBP1-85277

Dynactin Subunit 2/DCTN2/DCTN-50 Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Host	Rabbit
Gene ID	10540
Gene Symbol	DCTN2
Species	Human, Mouse, Rat
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: RWSPIASTLPELVQRLVTIKQLHEQAMQFGQLLTHLDTTQQMIANSLKDNTLLT QVQTTMRENLATVEGNFASIDERMKKL


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Knockdown Validated
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:2500 - 1:5000, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:2500 - 1:5000, Knockdown Validated
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF and WB reactivity reported in (PMID: 25669897).

Images


Western Blot: Dynactin Subunit 2/DCTN2/DCTN-50 Antibody [NBP1-85277] - Quantification of Dynactin-2 knockdown in HeLa cells by Western Blot. Western blots were cropped as indicated (black boxes). For all samples, SDS-PAGE and Western blotting were performed under the same experimental conditions. n:20 for each condition. * = $p < 0.05$; ** = $p < 0.01$. Scale bar 10 um. Error bars represent SEM. Image collected and cropped by CiteAb from the following publication (<https://www.nature.com/articles/srep08389>) licensed under a CC-BY license.




Immunohistochemistry-Paraffin: Dynactin Subunit 2/DCTN2/DCTN-50 Antibody [NBP1-85277] - Staining of human cerebral cortex, colon, kidney and testis using Anti-DCTN2 antibody NBP1-85277 (A) shows similar protein distribution across tissues to independent antibody NBP1-85278 (B).




Western Blot: Dynactin Subunit 2/DCTN2/DCTN-50 Antibody [NBP1-85277] - Analysis using Anti-DCTN2 antibody NBP1-85277 (A) shows similar pattern to independent antibody NBP1-85278 (B).




Immunohistochemistry-Paraffin: Dynactin Subunit 2/DCTN2/DCTN-50 Antibody [NBP1-85277] - Staining of human kidney.




Western Blot: Dynactin Subunit 2/DCTN2/DCTN-50 Antibody [NBP1-85277] - Lane 1: NIH-3T3 cell lysate (Mouse embryonic fibroblast cells). Lane 2: NBT-II cell lysate (Rat Wistar bladder tumor cells).




Immunohistochemistry-Paraffin: Dynactin Subunit 2/DCTN2/DCTN-50 Antibody [NBP1-85277] - Staining of human testis.



Immunohistochemistry-Paraffin: Dynactin Subunit 2/DCTN2/DCTN-50 Antibody [NBP1-85277] - Staining of human cerebral cortex.



Immunohistochemistry-Paraffin: Dynactin Subunit 2/DCTN2/DCTN-50 Antibody [NBP1-85277] - Staining of human colon.



Publications

Morris EJ, Assi K, Salh B, Dedhar S. Integrin-Linked Kinase links Dynactin-1/Dynactin-2 with cortical Integrin receptors to orient the mitotic spindle relative to the substratum Sci Rep. 2015-02-11 [PMID: 25669897] (ICC/IF, WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-85277

NBP1-85277PEP	Dynactin Subunit 2/DCTN2/DCTN-50 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-85277

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

