

# Product Datasheet

## Rab13 Antibody - BSA Free

### NBP1-85799

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

#### Publications: 5

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-85799](http://www.novusbio.com/NBP1-85799)

Updated 3/4/2026 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-85799](http://www.novusbio.com/reviews/destination/NBP1-85799)



**NBP1-85799**

Rab13 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Host	Rabbit
Gene ID	5872
Gene Symbol	RAB13
Species	Human, Mouse
Reactivity Notes	Use in Mouse reported in scientific literature (PMID:31290739).
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: AMGIILVYDITDEKSFENIQNWKSIKENASAGVERLLLGNKCDMEAKRKVQKE QADKLAREHGIRFFETSAKSSMNVDFAFSSLARDILLKSGGRRSGNGNKPPST DLKTCDKKNTNKCS


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Proximity Ligation Assay, Knockdown Validated
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:50 - 1:200, Proximity Ligation Assay Reported in scientific literature (PMID:31290739), Knockdown Validated
Application Notes	IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF, Fixation Permeabilization: Use PFA/Triton X-100.

**Images**

Rab13-Antibody-Knockdown-Validated-NBP1-85799-img0015.jpg




Rab13-Antibody-Knockdown-Validated-NBP1-85799-img0016.jpg




Immunocytochemistry/Immunofluorescence: Rab13 Antibody [NBP1-85799] - Staining of human cell line U-2 OS shows localization to plasma membrane and cytosol. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: Rab13 Antibody [NBP1-85799] - Staining of human testis shows moderate membranous positivity in cells in seminiferous ducts.




Immunohistochemistry-Paraffin: Rab13 Antibody [NBP1-85799] - Staining of human stomach shows moderate membranous positivity in glandular cells.




Immunohistochemistry-Paraffin: Rab13 Antibody [NBP1-85799] - Staining of human prostate shows weak membranous positivity in glandular cells.




Immunohistochemistry-Paraffin: Rab13 Antibody [NBP1-85799] - Staining of human skeletal muscle shows no positivity in myocytes as expected.




Western Blot: Rab13 Antibody [NBP1-85799] - RABIF distribution & effect on RAB13 expression  
 Representative immunofluorescence images of RABIF protein in cells transfected with the indicated siRNAs. Reduction of intensity in knockdown cells confirms the specificity of the signal. Calibration bar shows intensity values. Note that RABIF exhibits a mostly perinuclear enrichment. RAB13 expression in cells with CRISPR knockdown of RABIF using the indicated sgRNAs (see also Fig 7C). Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/32946136>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.




Western Blot: Rab13 Antibody [NBP1-85799] - RAB13 RNA mislocalization does not affect RAB13 binding to membranes or association with REP1 or RabGDI  
 Cells treated with the indicated PMOs were fractionated into soluble & particulate fractions, & the indicated proteins were detected by Western blot. RhoGDI & TfRc serve as soluble & particulate markers, respectively. Lysates from the indicated GFP or GFP-RAB13-expressing cell lines were immunoprecipitated with anti-GFP antibodies & blotted to detect the indicated proteins. Relative REP1 & RabGDI binding are quantified in the graphs from  $n = 3$  (REP1) &  $n = 5$  (RabGDI) independent experiments. Bars: mean  $\pm$  s.e.m. No significant differences were detected by Wilcoxon matched-pairs signed-rank test. Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/32946136>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.




Western Blot: Rab13 Antibody [NBP1-85799] - RAB13 RNA & protein exhibit distinct subcellular distributions. Representative FISH images showing RAB13 RNA distribution in MDA-MB-231 cells. Nuclei & cell outlines are shown in blue & green, respectively. Arrows point to RAB13 RNA concentrated at protrusive regions. Boxed regions are magnified in the insets. Representative immunofluorescence images of RAB13 protein in cells transfected with the indicated siRNAs. Reduction of intensity in RAB13 knockdown cells confirms the specificity of the signal. Arrows point to perinuclear RAB13 protein. Calibration bar shows intensity values. Ratios of peripheral/perinuclear intensity calculated from images as shown in (A) & (B). Bars: mean  $\pm$  s.e.m. Values within each bar represent number of cells observed in 3 independent experiments. Protrusions (Ps) & cell bodies (CB) of cells induced to migrate toward LPA were isolated & analyzed to detect the indicated proteins (by Western blot; left panels) or RNAs (by RT-ddPCR; right panel). Ps/CB enrichment ratios from 2 independent experiments are shown. Bars: mean  $\pm$  s.e.m. The enrichment of pY397-FAK serves to verify the enrichment of protrusions containing newly formed adhesions in the Ps fraction. Data information: P-values: \*\* < 0.01; \*\*\*\* < 0.0001 by Student's t-test (C) or analysis of variance with Dunnett's multiple comparisons test against GFP (D). Scale bars: 10  $\mu$ m. Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/32946136>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Western Blot: Rab13 Antibody [NBP1-85799] - Peripheral RAB13 RNA translation is required for RAB13 protein activation but not steady-state distribution or membrane association. Wide-field images of RAB13 immunofluorescence in MDA-MB-231 cells treated with control or RAB13 (191 + 230) PMOs & ratios of peripheral/perinuclear intensity. Scale bars: 10  $\mu$ m. n = 45–50 cells. Bars: mean  $\pm$  s.e.m. Similar results were observed in two additional independent experiments. Fluorescence images (projections of confocal slices throughout the cell height) of cells expressing GFP-RAB13 with the indicated UTRs. Note that in both cases the protein assumes indistinguishable distribution. Scale bars: 10  $\mu$ m. Soluble/particulate fractionation of the indicated cell lines followed by Western blot to detect the indicated proteins. RhoGDI & TfRc serve as soluble & particulate markers, respectively. Graph shows quantitation from n = 3 independent experiments. Bars: mean  $\pm$  s.e.m. D, E Active RAB13 (RAB13-GTP) was pulled down using MICAL-L1 RBD-GST from the indicated PMO-treated cells (D) or GFP-RAB13-expressing lines (E). The amount of endogenous or exogenous RAB13 was measured by quantitative Western blot, & relative levels of active RAB13 are plotted. n = 8 (D), n = 4 (E). Bars: mean  $\pm$  s.e.m. Data information: P-values: \* < 0.05, \*\* < 0.01 by Kruskal-Wallis test. Source data are available online for this figure. Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/32946136>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.




**Immunocytochemistry/ Immunofluorescence: Rab13 Antibody [NBP1-85799] - RAB13 protein levels do not change upon serum stimulation** MDA-MB-231 cells were stimulated with serum for the indicated times, in the presence or absence of cycloheximide (CHX). The cells were also treated with control or RAB13 mislocalizing PMOs. Representative Western blot analysis of whole cell lysates & corresponding quantitations of RAB13 levels from n = 4–5 replicates. Bars: mean  $\pm$  s.e.m. No significant differences by Friedman's test. Increase in pY397-FAK levels attests to serum stimulation. RAB13 immunofluorescence at representative protrusive regions. A cell line expressing GFP was used to delineate cell borders & provide an internal cytosolic control. RAB13 signal at front lamellipodial regions was quantified. n = 35–51 protrusions. Bars: mean  $\pm$  s.e.m. No increase is detected upon stimulation. By contrast, at early time points a decrease is detected (5 & 20 min,  $P < 0.01$  by Kruskal–Wallis test), potentially arising from serum-induced endocytosis of RAB13-containing membranes. Scale bars: 8  $\mu$ m. Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/32946136>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



**Western Blot: Rab13 Antibody [NBP1-85799] - RAB13 protein levels do not change upon serum stimulation** MDA-MB-231 cells were stimulated with serum for the indicated times, in the presence or absence of cycloheximide (CHX). The cells were also treated with control or RAB13 mislocalizing PMOs. Representative Western blot analysis of whole cell lysates & corresponding quantitations of RAB13 levels from n = 4–5 replicates. Bars: mean  $\pm$  s.e.m. No significant differences by Friedman's test. Increase in pY397-FAK levels attests to serum stimulation. RAB13 immunofluorescence at representative protrusive regions. A cell line expressing GFP was used to delineate cell borders & provide an internal cytosolic control. RAB13 signal at front lamellipodial regions was quantified. n = 35–51 protrusions. Bars: mean  $\pm$  s.e.m. No increase is detected upon stimulation. By contrast, at early time points a decrease is detected (5 & 20 min,  $P < 0.01$  by Kruskal–Wallis test), potentially arising from serum-induced endocytosis of RAB13-containing membranes. Scale bars: 8  $\mu$ m. Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/32946136>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



**Analysis in U-87MG ATCC cells transfected with control siRNA, target specific siRNA probe #1 and #2, using Anti-RAB13 antibody. Remaining relative intensity is presented. Loading control: Anti-GAPDH.**



**Publications**

Moissoglu K, Stueland M, Gasparski AN et al. RNA localization and co-translational interactions control RAB13 GTPase function and cell migration EMBO J. 2020-09-18 [PMID: 32946136] (ICC/IF, WB, Human)

Translational regulation of protrusion-localized RNAs involves silencing and clustering after transport Moissoglu K, Yasuda K, Wang T et al. Elife [PMID: 31290739] (PLA, WB, Mouse)

Eaton A Effects of Stretch on the Bladder Umbrella Cell Apical Junctional Ring Thesis 2020-01-01 (WB, Rat)

Eaton AF, Clayton DR, Ruiz WG et al. Expansion and Contraction of the Umbrella Cell Apical Junctional Ring in Response to Bladder Filling and Voiding Mol. Biol. Cell 2019-06-05 [PMID: 31166831]

Abou-Zeid N, Pandjaitan R, Sengmanivong L et al. MICAL-like1 mediates epidermal growth factor receptor endocytosis. Mol Biol Cell 2011-09-01 [PMID: 21795389]





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-85799**

---

NBP1-85799PEP	Rab13 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-85799](http://www.novusbio.com/reviews/submit/NBP1-85799)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

