

Product Datasheet

PYGL Antibody - BSA Free

NBP1-86147

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-86147

Updated 12/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-86147



NBP1-86147

PYGL Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Description	Novus Biologicals Rabbit PYGL Antibody - BSA Free (NBP1-86147) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5836
Gene Symbol	PYGL
Species	Human
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: VVAATLQDIIRRFKASKFGSTRGAGTVFDAFPDQVAIQLNDTHPALAIPELMRIFV DIEKLPWSKAWELTQKTFAYTNHTVLPEALERWPVDLVEKLLPRHLEIYEINQK HLDRIVALFPKDVDRLRMSLI


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:200 - 1:500, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:200 - 1:500
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100.

Images


Analysis in human cell lines PC-3 and SK-MEL-30 using Anti-PYGL antibody. Corresponding PYGL RNA-seq data are presented for the same cell lines. Loading control: Anti-PFN1.




Immunohistochemistry-Paraffin: PYGL Antibody [NBP1-86147] - Analysis in human liver and pancreas tissues. Corresponding PYGL RNA-seq data are presented for the same tissues.




Immunohistochemistry-Paraffin: PYGL Antibody [NBP1-86147] - Staining of human liver shows strong cytoplasmic positivity in hepatocytes.




Immunohistochemistry-Paraffin: PYGL Antibody [NBP1-86147] - Staining of human pancreas shows very weak positivity in exocrine glandular cells.




Immunohistochemistry-Paraffin: PYGL Antibody [NBP1-86147] - Staining of human bone marrow shows moderate cytoplasmic positivity in hematopoietic cells.




Immunohistochemistry-Paraffin: PYGL Antibody [NBP1-86147] - Staining of human cerebral cortex shows very weak positivity in neuronal cells.




Immunohistochemistry-Paraffin: PYGL Antibody [NBP1-86147] - Staining of human testis shows very weak positivity in cells in seminiferous ducts.



Immunohistochemistry-Paraffin: PYGL Antibody [NBP1-86147] - Staining of human bone marrow, cerebral cortex, liver and testis using Anti-PYGL antibody (A) NBP1-86147 shows similar protein distribution across tissues to independent antibody NBP1-86182 (B).



Staining of human cell line U-2 OS shows localization to plasma membrane & cytosol.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-86147

NBP1-86147PEP	PYGL Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-86147

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

