

Product Datasheet

NCAPG2 Antibody - BSA Free

NBP1-86633

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-86633

Updated 12/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-86633



NBP1-86633

NCAPG2 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol
Target Molecular Weight	131 kDa

Product Description	
Description	Novus Biologicals Rabbit NCAPG2 Antibody - BSA Free (NBP1-86633) is a polyclonal antibody validated for use in IHC and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	54892
Gene Symbol	NCAPG2
Species	Human
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: QVGHILELVDNWLPTTEHAQAKSNTASKGRVQIHDTRPVKPELALVYIEYLLTHP KNRECLLSAPRKLNHLLKALETSKADLESLLQTPGGKPRGFSEAAAPRAFGLH CRLSIHLQH

Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:50 - 1:200
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF, Fixation/Permeabilization: PFA/Triton X-100.

Images


Immunohistochemistry-Paraffin: NCAPG2 Antibody [NBP1-86633] - Staining of human tonsil shows moderate nuclear positivity in germinal center and non-germinal center cells.




Immunohistochemistry-Paraffin: NCAPG2 Antibody [NBP1-86633] - Staining of human testis shows moderate to strong nuclear positivity in cells in seminiferous ducts and Leydig cells.




Immunohistochemistry-Paraffin: NCAPG2 Antibody [NBP1-86633] - Staining of human prostate shows moderate nuclear positivity in glandular cells.



Immunohistochemistry-Paraffin: NCAPG2 Antibody [NBP1-86633] - Staining of human skeletal muscle shows weak to moderate nuclear positivity in myocytes.



Staining of human cell line A-431 shows localization to nuclear speckles.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-86633

NBP1-86633PEP	NCAPG2 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-86633

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

