

Product Datasheet

BRD4 Antibody - BSA Free

NBP1-86640

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-86640

Updated 2/1/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-86640



NBP1-86640

BRD4 Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol


Product Description	
Description	Novus Biologicals Rabbit BRD4 Antibody - BSA Free (NBP1-86640) is a polyclonal antibody validated for use in IHC, ICC/IF, Simple Western and ChIP. Anti-BRD4 Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	23476
Gene Symbol	BRD4
Species	Human
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: PQPAKPQQVIQHSHSPRHHKSDPYSTGHLREAPSPLMIHSPQMSQFQSLTHQ SPPQQNVQPKKQELRAASVVQPQLVVKKEEKIHSPIIRSEPFSPSLRPEPPKH PESIKAPVHLPQRPEMKPVDVGRPVI RPPEQNAPP

Product Application Details	
Applications	Simple Western, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Chromatin Immunoprecipitation (ChIP)
Recommended Dilutions	Simple Western 1:60, Immunohistochemistry 1:500 - 1:1000, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:500 - 1:1000, Chromatin Immunoprecipitation (ChIP) Reported in scientific literature (PMID:25616107).
Application Notes	ICC/IF, Fixation Permeabilization: Use PFA/Triton X-100. IHC-Paraffin, HIER pH6 retrieval is recommended. Simple Western only 10 - 15 uL of the recommended dilution is used per data point. See Simple Western Antibody Database for Simple Western validation: Tested in U-251MG sp, separated by Size, antibody dilution of 1:60, apparent MW was 274 kDa. Separated by Size-Wes, Sally Sue/Peggy Sue.




Images


Staining of human cell line U-251 MG shows localization to nucleoplasm.




Immunohistochemistry-Paraffin: BRD4 Antibody [NBP1-86640] - Staining of human tonsil shows moderate to strong nuclear positivity in non-germinal center cells.




Immunohistochemistry-Paraffin: BRD4 Antibody [NBP1-86640] - Staining of human cerebral cortex shows moderate to strong nuclear positivity in neurons.



Immunohistochemistry-Paraffin: BRD4 Antibody [NBP1-86640] - Staining of human placenta shows moderate to strong nuclear positivity in trophoblastic cells.




Immunohistochemistry-Paraffin: BRD4 Antibody [NBP1-86640] - Staining of human testis shows moderate to strong nuclear positivity in cells in seminiferous ducts.



Simple Western: BRD4 Antibody [NBP1-86640] - Simple Western lane view shows a specific band for BRD4 in 0.2 mg/ml of U-251MG sp lysate. This experiment was performed under reducing conditions using the 66-440 kDa separation system.



Simple Western: BRD4 Antibody [NBP1-86640] - Electropherogram image(s) of corresponding Simple Western lane view. BRD4 antibody was used at 1:60 dilution on U-251MG sp lysates(s).



Publications

Zimmermann M, Aguilera FB, Castellucci M et al. Chromatin remodelling and autocrine TNFalpha are required for optimal interleukin-6 expression in activated human neutrophils. *Nat Commun.* [PMID: 25616107] (Chemotaxis, Human)

Details:

Citation using the PE form of this antibody.

Zuber J, Shi J, Wang E et al. RNAi screen identifies Brd4 as a therapeutic target in acute myeloid leukaemia. *Nature* 2011-08-01 [PMID: 21814200]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-86640

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-38931PEP	BRD4 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-86640

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

