

Product Datasheet

RABL4 Antibody - BSA Free NBP1-87170

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 4

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-87170

Updated 12/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-87170



NBP1-87170

RABL4 Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol


Product Description	
Description	Novus Biologicals Rabbit RABL4 Antibody - BSA Free (NBP1-87170) is a polyclonal antibody validated for use in WB. Anti-RABL4 Antibody: Cited in 4 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	11020
Gene Symbol	IFT27
Species	Human
Reactivity Notes	Rat (80%).
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: TNEESFNCSKWLEKARSQAPGISLPGVLVGNKTDLAGRRAVDSAEARAWAL GQGLECFETSVKEMENFEAPFHCLAKQFHQLYREKVEV

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml




Images


Western Blot: RABL4 Antibody [NBP1-87170] - Analysis in human cell lines HEK293 and Caco-2 using anti-IFT27 antibody. Corresponding IFT27 RNA-seq data are presented for the same cell lines. Loading control: anti-COX4I1.




Western Blot: RABL4 Antibody [NBP1-87170] - IFT proteins associate with central spindle MTs in anaphase and are required for their organization. Western blots (left) showing the amount of IFT27 and IFT88 in LLC-PK1 cells transfected with control (CT), IFT27 or IFT88 siRNA. alpha-tubulin: loading control. Immunofluorescence images (right) showing a decrease of IFT88 staining in an anaphase cell upon IFT88 depletion. Scale bars: 10?um for main images, 3?um for insets. Image collected and cropped by CiteAb from the following publication (<https://www.nature.com/articles/s41467-017-01479-3>), licensed under a CC-BY license.



Western Blot: RABL4 Antibody [NBP1-87170] - Analysis in control (vector only transfected HEK293T lysate) and IFT27 over-expression lysate (Co-expressed with a C-terminal myc-DDK tag (3.1 kDa) in mammalian HEK293T cells).



Western Blot: RABL4 Antibody [NBP1-87170] - IFT proteins are required for the timely relocalization of the kinase Aurora B to the central spindle. Myc immunoprecipitation performed on mitotic HeLa Kyoto cells transfected with Myc or with Myc-Aurora B. Image collected and cropped by CiteAb from the following publication (<https://www.nature.com/articles/s41467-017-01479-3>), licensed under a CC-BY license.



Publications

kanie T, Abbott kL, MooneY NA et al. The CEP19-RABL2 GTPase Complex Binds IFT-B to Initiate Intraflagellar Transport at the Ciliary Base. *Dev. Cell.* 2017-07-10 [PMID: 28625565]

Vitre B, Taulet N et al. IFT proteins interact with HSET to promote supernumerary centrosome clustering in mitosis. *EMBO Rep* 2020-04-06 [PMID: 32270908] (WB, Human)

Taulet N, Vitre B, Anguille C et al. IFT proteins spatially control the geometry of cleavage furrow ingression and lumen positioning. *Nat Commun* 2017-12-04 [PMID: 29203870] (WB)

Stadler C, Rexhepaj E, Singan VR et al. Immunofluorescence and fluorescent-protein tagging show high correlation for protein localization in mammalian cells. *Nat Methods* 2013-04-01 [PMID: 23435261]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-87170

NBP1-87170PEP	RABL4 Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-87170

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

