

Product Datasheet

FAM118A Antibody - BSA Free NBP1-88590

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-88590

Updated 4/8/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-88590



NBP1-88590

FAM118A Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Host	Rabbit
Gene ID	55007
Gene Symbol	FAM118A
Species	Human
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: VTQDAEVMEVLQNLRYRTKSFLVGCGETLRDQIFQALFLYSVPNKVDLEHYML VLKENEDHFFKHQADMLLHGKVVSYGDCFDHFPGYVQDLATQICKQQSPDAD RVDSTLLGNACQDCAKRKLEENGIE

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:50-1:200
Application Notes	IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF, Fixation Permeabilization: Use PFA/Triton X-100.

Images


Staining of human cell line U-2 OS shows localization to nucleus, cytosol & intermediate filaments.




Immunohistochemistry-Paraffin: FAM118A Antibody [NBP1-88590] - Staining of human endometrium shows strong cytoplasmic positivity in glandular cells.




Immunohistochemistry-Paraffin: FAM118A Antibody [NBP1-88590] - Staining of human lung shows strong cytoplasmic positivity in macrophages.




Immunohistochemistry-Paraffin: FAM118A Antibody [NBP1-88590] - Staining of human cerebellum shows strong cytoplasmic positivity in Purkinje cells.



Staining of human liver shows no positivity in hepatocytes as expected.



Analysis in control (vector only transfected HEK293T lysate) and FAM118A over-expression lysate (Co-expressed with a C-terminal myc-DDK tag (~3.1 kDa) in mammalian HEK293T cells).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-88590

NBP1-88590PEP	FAM118A Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-88590

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

