

Product Datasheet

PRKAG2 Antibody - BSA Free

NBP1-89324

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 4

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-89324

Updated 12/2/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-89324



NBP1-89324

PRKAG2 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Description	Novus Biologicals Rabbit PRKAG2 Antibody - BSA Free (NBP1-89324) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-PRKAG2 Antibody: Cited in 3 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	51422
Gene Symbol	PRKAG2
Species	Human, Mouse, Rat
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: STPTQVTKQHTFPLESYKHEPERLENRIYASSSPDGTGQRFCPSSFQSPTRPP LASPTHYAPSKAAALAAALGPAEAGMLEKLEFEDEAVEDSESGVYMRFRSHK


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:50 - 1:200
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100.

Images


Western Blot: PRKAG2 Antibody [NBP1-89324] - Western blot analysis in mouse cell line NIH-3T3, rat cell line NBT-II and rat cell line pC12.




Immunocytochemistry/Immunofluorescence: PRKAG2 Antibody [NBP1-89324] - Staining of human cell line U-2 OS shows localization to nucleoplasm. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: PRKAG2 Antibody [NBP1-89324] - Staining of human prostate shows moderate to strong positivity in nuclear membrane in glandular cells.




Western Blot: PRKAG2 Antibody [NBP1-89324] - Analysis in human cell line A-431.




Immunohistochemistry-Paraffin: PRKAG2 Antibody [NBP1-89324] - Staining of human skin shows moderate to strong positivity in nuclear membrane in keratinocytes



Immunohistochemistry-Paraffin: PRKAG2 Antibody [NBP1-89324] - Staining of human kidney shows moderate to strong positivity in nuclear membrane in cells in tubules.



Immunohistochemistry-Paraffin: PRKAG2 Antibody [NBP1-89324] - Staining of human testis shows moderate to strong positivity in nuclear membrane in cells in seminiferous ducts.



Publications

Li M, Wei X, Xiong J et al. Hierarchical inhibition of mTORC1 by glucose starvation-triggered AXIN lysosomal translocation and by AMPK Life Metabolism 2023-03-01 (Western Blot)

Tobias IS, Lazauskas KK, Arevalo JA et al. Fiber type-specific analysis of AMPK isoforms in human skeletal muscle: advancement in methods via capillary nano-immunoassay. Journal of Applied Physiology 2017-12-14 [PMID: 29357518] (Human)

Pinter K, Jefferson A, Czibik G et al. Subunit composition of AMPK trimers present in the cytoskeletal apparatus: Implications for drug target identification. Cell Cycle 2012-03-01 [PMID: 22333580]

Martinez T, Frances S, Lopez JM et al. Generation of digital responses in stress sensors. J Biol Chem 2009-09-01 [PMID: 19570986]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-89324

NBP1-89324PEP	PRKAG2 Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-89324

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

