

Product Datasheet

Human P-Selectin/CD62P ELISA Kit (Colorimetric) NBP1-89873

Unit Size: 96 Tests

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-89873

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-89873



NBP1-89873**Human P-Selectin/CD62P ELISA Kit (Colorimetric)****Product Information**

Unit Size	96 Tests
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.

Product Description


Description	This kit is an enzyme-linked immunosorbent assay for quantitative detection.
Gene ID	6403
Gene Symbol	SELP
Species	Human
Specificity/Sensitivity	The cross reactivity and interference of circulating factors of the immune system was evaluated by spiking these proteins at physiologically relevant concentrations into a human P-selectin positive sample. There was no cross reactivity or interference detected, notably not with human E-selectin and human L-selectin.
Kit Components	Aluminium pouch(es) with a Microwell Plate coated with monoclonal antibody to GLP-1
Standard Curve Range	0.63 - 40.00 ng/mL
Sensitivity	0.20 ng/mL
Inter Assay	CV% < 5.4%
Intra Assay	CV% < 7.8%
Assay Type	Sandwich ELISA
Spiking Recovery	For recovery data see the manual.
Suitable Sample Type	Cell culture supernatant, serum, plasma (EDTA, citrate, heparin), amniotic fluid
Sample Volume	10 uL

Product Application Details

Applications	ELISA
Recommended Dilutions	ELISA

Images

ELISA: Human P-Selectin/CD62P ELISA Kit (Colorimetric) [NBP1-89873] - Representative standard curve for human P-selectin ELISA. Human P-selectin was diluted in serial 2-fold steps in Sample Diluent. Do not use this standard curve to derive test results. A standard curve must be run for each group of microwell strips assayed.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-89873

NBP1-85744PEP	P-Selectin/CD62P Recombinant Protein Antigen
210-TA-005	TNF-alpha [Unconjugated]
137-PS-050	P-Selectin/CD62P
M6000B-1	IL-6 [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-89873

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

