

Product Datasheet

Fetuin A/AHSG Antibody - BSA Free NBP1-90302

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-90302

Updated 12/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-90302



NBP1-90302


Fetuin A/AHSG Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol
Product Description	
Description	Novus Biologicals Rabbit Fetuin A/AHSG Antibody - BSA Free (NBP1-90302) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-Fetuin A/AHSG Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	197
Gene Symbol	AHSG
Species	Human, Porcine
Reactivity Notes	Porcine reactivity reported in scientific literature (PMID: 24862986).
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: GLIYRQPNCDDEPETEEAALVAIDYINQNLPGWKHTLNQIDEVKVWPQQPSGE LFEIEIDTLETTCHVLDPTPVARCSVRQLKEHAVEGDCDFQLLKLDGKFSVVYAK CDSSPDSAEDVRKVCQDCPLLAPLNDTRVVHAAKAALA
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Knockdown Validated
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:50 - 1:200, Knockdown Validated
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100.




Images

Staining of human cell line Hep G2 shows localization to the Golgi apparatus.




Immunohistochemistry-Paraffin: Fetuin A/AHSG Antibody [NBP1-90302]
 - Staining of human bone marrow, cerebral cortex, liver and lymph node using Anti-AHSG antibody NBP1-90302 (A) shows similar protein distribution across tissues to independent antibody NBP1-90303 (B).




Immunohistochemistry-Paraffin: Fetuin A/AHSG Antibody [NBP1-90302]
 - Staining of human bone marrow.




Immunohistochemistry-Paraffin: Fetuin A/AHSG Antibody [NBP1-90302]
 - Staining of human liver.




Immunohistochemistry-Paraffin: Fetuin A/AHSG Antibody [NBP1-90302]
- Staining of human cerebellum shows moderate cytoplasmic positivity in a subset of Purkinje cells.




Immunohistochemistry-Paraffin: Fetuin A/AHSG Antibody [NBP1-90302]
- Staining of human cerebral cortex.




Immunohistochemistry-Paraffin: Fetuin A/AHSG Antibody [NBP1-90302]
- Staining of human kidney shows moderate extra-cellular positivity and moderate cytoplasmic positivity in cells in tubules.




Staining of human skeletal muscle shows no cytoplasmic positivity in myocytes as expected.




Immunohistochemistry-Paraffin: Fetuin A/AHSG Antibody [NBP1-90302]
- Staining of human lymph node.



Staining of human stomach shows moderate positivity in extracellular matrix.



Analysis in Caco-2 cells transfected with control siRNA, target specific siRNA probe #1 and #2. Remaining relative intensity is presented. Loading control: Anti-GAPDH.



Publications

Hossain MA, De Souza AI, Bagul A et al. HSP70, Peroxiredoxin-3 and -6 are upregulated during renal warm ischaemia in a donation after circulatory death model. J Proteomics. 2014-05-24 [PMID: 24862986] (WB, Porcine)

Levin Y, Wang L, Schwarz E et al. Global proteomic profiling reveals altered proteomic signature in schizophrenia serum. Mol Psychiatry 2010-11-01 [PMID: 19546861]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-90302

NBP1-90302PEP	Fetuin A/AHSG Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-90302

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

