

# Product Datasheet

## Fibrinogen beta chain Antibody - BSA Free NBP1-90955

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

### Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-90955](http://www.novusbio.com/NBP1-90955)

Updated 12/2/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-90955](http://www.novusbio.com/reviews/destination/NBP1-90955)



**NBP1-90955**

Fibrinogen beta chain Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol


Product Description	
Description	Novus Biologicals Rabbit Fibrinogen beta chain Antibody - BSA Free (NBP1-90955) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-Fibrinogen beta chain Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	2244
Gene Symbol	FGB
Species	Human
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: VIQNRQDGSVDFGRKWDPYKQGFGNVATNTDGKNYCGLPGEYWLGNKISQ LTRMGPTELLIEMEDWKGDVKVKAHYGGFTVQNEANKYQISVNKYRGTAGNAL MDGASQLMGENRTMTIHNGMFFSTYDRDNDGWLTSDPRKQC

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:1000 - 1:2500, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:1000 - 1:2500
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF, Fixation/Permeabilization: PFA/Triton X-100.




## Images


Immunohistochemistry-Paraffin: Fibrinogen beta chain Antibody [NBP1-90955] - Staining of human cerebral cortex, liver, placenta and rectum using Fibrinogen beta chain antibody NBP1-90955 (A) shows similar protein distribution across tissues to independent antibody NBP1-90956 (B).




Western Blot: Fibrinogen beta chain Antibody [NBP1-90955] - Analysis using Fibrinogen beta chain antibody NBP1-90955 (A) shows similar pattern to independent antibody NBP1-90956 (B).




Immunocytochemistry/Immunofluorescence: Fibrinogen beta chain Antibody [NBP1-90955] - Staining of human cell line Hep G2 shows localization to endoplasmic reticulum.




Immunohistochemistry-Paraffin: Fibrinogen beta chain Antibody [NBP1-90955] - Staining of human liver shows weak cytoplasmic positivity in hepatocytes.



Immunohistochemistry-Paraffin: Fibrinogen beta chain Antibody [NBP1-90955] - Staining of human cerebral cortex shows no positivity in neurons as expected.



Immunohistochemistry-Paraffin: Fibrinogen beta chain Antibody [NBP1-90955] - Staining of human placenta shows strong positivity in plasma in blood vessels.



Immunohistochemistry-Paraffin: Fibrinogen beta chain Antibody [NBP1-90955] - Staining of human rectum shows moderate positivity in plasma in blood vessels, as well as positivity in extracellular matrix.



## Publications

Qundos U, Hong MG, Tybring G et al. Profiling post-centrifugation delay of serum and plasma with antibody bead arrays. J Proteomics 2013-12-16 [PMID: 23631827]

Kato BS, Nicholson G, Neiman M et al. Variance decomposition of protein profiles from antibody arrays using a longitudinal twin model. Proteome Sci 2011-11-17 [PMID: 22093360]



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-90955**

---

NBP1-90955PEP	Fibrinogen beta chain Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-90955](http://www.novusbio.com/reviews/submit/NBP1-90955)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

