

Product Datasheet

PKC delta Antibody - BSA Free NBP1-90957

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-90957

Updated 12/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-90957



NBP1-90957

PKC delta Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Description	Novus Biologicals Rabbit PKC delta Antibody - BSA Free (NBP1-90957) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5580
Gene Symbol	PRKCD
Species	Human
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: AFNSYELGSLQAEDEANQPFCVAKMKEALSTERGKTLVQKKPTMYPEWKSTF DAHIYEGRVIQIVLMRAAEVPSEVTGVSVLAERCKKNNNGKAEFWLDLQPQA KVLMSVQYFLEDVDCKQSMRSEDEAKFPTMNRRAIKQ


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:20 - 1:50, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:20-1:50
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100.

Images


Western Blot: PKC delta Antibody [NBP1-90957] - Lane 1: Marker [kDa] 220, 112, 84, 47, 32, 26, 17. Lane 2: Human cell line RT-4



Immunocytochemistry/Immunofluorescence: PKC delta Antibody [NBP1-90957] - Staining of human cell line A-431 shows localization to cytosol. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: PKC delta Antibody [NBP1-90957] - Staining of human upper gastrointestinal shows moderate cytoplasmic positivity in glandular cells.




Immunohistochemistry-Paraffin: PKC delta Antibody [NBP1-90957] - Staining of human prostate shows high expression.




Immunohistochemistry-Paraffin: PKC delta Antibody [NBP1-90957] - Staining of human skeletal muscle shows low expression as expected.




Immunohistochemistry-Paraffin: PKC delta Antibody [NBP1-90957] - Staining in human prostate and skeletal muscle tissues using anti-PRKCD antibody. Corresponding PRKCD RNA-seq data are presented for the same tissues.




Immunohistochemistry-Paraffin: PKC delta Antibody [NBP1-90957] - Analysis in human tonsil and skeletal muscle tissues. Corresponding PKC delta RNA-seq data are presented for the same tissues.




Immunohistochemistry-Paraffin: PKC delta Antibody [NBP1-90957] - Staining of human prostate, skeletal muscle, tonsil and upper gastrointestinal using Anti-PKC delta antibody NBP1-90957 (A) shows similar protein distribution across tissues to independent antibody NBP1-90351 (B).




Immunohistochemistry-Paraffin: PKC delta Antibody [NBP1-90957] - Staining of human prostate shows moderate cytoplasmic positivity in glandular cells.



Immunohistochemistry-Paraffin: PKC delta Antibody [NBP1-90957] - Staining of human skeletal muscle shows only very weak cytoplasmic positivity in myocytes.




Immunohistochemistry: PKC delta Antibody [NBP1-90957] - Staining of human tonsil shows strong cytoplasmic positivity in germinal center cells.



Lane 1: Marker [kDa] 220, 112, 84, 47, 32, 26, 17

Lane 2: Human cell line RT-4





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-90957

NBP1-90957PEP	PKC delta Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-90957

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

