

Product Datasheet

Schlafen 11 Antibody - BSA Free NBP1-92368

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 3

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-92368

Updated 6/3/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-92368



NBP1-92368

Schlafen 11 Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Host	Rabbit
Gene ID	91607
Gene Symbol	SLFN11
Species	Human
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: REVLGCAKENVDPDSLRRKIEQAIYKLPCVHFCQPQRPITFTLKIVDVLKRGELY GYACMIRVNPFCFAVFSEAPNSWI

Product Application Details


Applications	Western Blot, Simple Western, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Simple Western 1:20, Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:50 - 1:200
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100. See Simple Western Antibody Database for Simple Western validation: Tested in MOLT-4, HEK 239T, separated by Size, antibody dilution of 1:20, apparent MW was 102 kDa

Images


Western Blot: Schlafen 11 Antibody [NBP1-92368] - Analysis in human cell lines HEK293 and Caco-2 using anti-SLFN11 antibody. Corresponding SLFN11 RNA-seq data are presented for the same cell lines. Loading control: anti-COX4I1.



Immunocytochemistry/Immunofluorescence: Schlafen 11 Antibody [NBP1-92368] - Staining of human cell line A-431 shows localization to nucleoplasm & cytosol. Antibody staining is shown in green.




Western Blot: Schlafen 11 Antibody [NBP1-92368] - Analysis in human cell line MOLT-4.



Simple Western: Schlafen 11 Antibody [NBP1-92368] - Simple Western lane view shows a specific band for SLFN11 in 0.2 mg/ml of MOLT-4 (left), HEK 239T (right) lysate. This experiment was performed under reducing conditions using the 12-230 kDa separation system.



Simple Western: Schlafen 11 Antibody [NBP1-92368] - Electropherogram image(s) of corresponding Simple Western lane view. Schlafen 8 antibody was used at 1:20 dilution on MOLT-4 and HEK 239T lysate(s).



Publications

Kim ET, Dybas JM, Kulej K, et al. Comparative proteomics identifies Schlafen 5 (SLFN5) as a herpes simplex virus restriction factor that suppresses viral transcription *Nature microbiology* 2021-02-01 [PMID: 33432153]

Stadler C, Rexhepaj E, Singan VR et al. Immunofluorescence and fluorescent-protein tagging show high correlation for protein localization in mammalian cells. *Nat Methods* 2013-04-01 [PMID: 23435261]

Li M, Kao E, Gao X et al. Codon-usage-based inhibition of HIV protein synthesis by human schlafen 11. *Nature* 2012-11-01 [PMID: 23000900]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-92368

NBP1-92368PEP	Schlafen 11 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-92368

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

