

Product Datasheet

EWSR1 Antibody (5H7) - BSA Free NBP1-92686

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 3


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-92686

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-92686



NBP1-92686

EWSR1 Antibody (5H7) - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	5H7
Preservative	5mM Sodium Azide
Isotype	IgG2b
Purity	Protein G purified
Buffer	50% PBS, 50% glycerol
Target Molecular Weight	80 kDa

Product Description	
Description	Novus Biologicals Mouse EWSR1 Antibody (5H7) - BSA Free (NBP1-92686) is a monoclonal antibody validated for use in IHC, WB and ICC/IF. Anti-EWSR1 Antibody: Cited in 3 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	2130
Gene Symbol	EWSR1
Species	Human, Mouse, Rat, Canine, Equine
Immunogen	Full length recombinant human EWSR1 expressed in and purified from E. coli. [UniProt# Q01844]


Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:1000-1:2000, Immunohistochemistry 1:1000, Immunocytochemistry/ Immunofluorescence 1:1000
Application Notes	This EWSR1 (5H7) antibody is useful for Immunocytochemistry/Immunofluorescence, Immunohistochemistry, and Western blot, where a band can be seen at approximately 80 kDa.

Images


Western Blot: EWSR1 Antibody (5H7) [NBP1-92686] - Image shows a specific band for EWSR1 (observed molecular weight ~95 kDa) in HEK293T and huh-7 lysate.




Immunocytochemistry/Immunofluorescence: EWSR1 Antibody (5H7) [NBP1-92686] - HeLa cell cultures stained with antibody NBP1-92686 (green) and chicken antibody to vimentin (NB300-223, red). Blue is a DNA stain, and it is clear that in these cells EWSR1 is localized along with the DNA in the nucleus.



Western Blot: EWSR1 Antibody (5H7) [NBP1-92686] - Analysis of different cell line lysates using EWSR1 antibody, dilution 1:1000 (Green): [1] protein standard (Red), [2] HeLa, [3] HEK293, [4] mouse NIH-3T3, [5] rat PC12, [6] equine NBL6 and [7] canine A72 cells. The strong band at ~80kDa corresponds to the EWS protein seen in all species tested.



Effects of siRNA-mediated knockdown of expression of AldoA, EWSR1 and ILF3-90 on enterovirus replication. A, B. HeLa cells were transfected with siRNAs specific to AldoA, EWSR1 and 90kDa isoform of ILF3, or non-targeting control siRNA, and polio or Coxsackie B3 replicon replication assays were performed 72 h post siRNA transfection. The total replication signal was calculated as the area under the corresponding kinetics curves. Cell viability signal is proportional to the level of ATP in cells. Western blots show the efficacy of siRNA-mediated knockdown of the targeted proteins. C. HeLa cells were transfected with siRNAs specific to AldoA, EWSR1 and 90kDa isoform of ILF3, or non-targeting control siRNA. 72 h post siRNA transfection cells were infected with an MOI of 1 PFU/cell of poliovirus or Coxsackie virus B3, and the total virus yield was determined at 6 h p. i. Western blots show the efficacy of siRNA-mediated knockdown of the targeted proteins. D. HeLa cells were transfected with siRNAs specific to 90kDa isoform of ILF3 or EWSR1 or a non-targeting control siRNA. 72h post siRNA transfection, the cells were transfected with a replication-defective replicon RNA containing the Δ 3D mutation. The total translation signal was calculated as the area under the curve from 1 to 4 h post Δ 3D RNA transfection. Western blots show the efficacy of siRNA-mediated knockdown of the targeted proteins. E. HeLa cells were transfected with siRNAs specific to AldoA, EWSR1 and 90kDa isoform of ILF3, or non-targeting control siRNA, and cell viability assays detecting the level of ATP or the activity of the mitochondrial respiratory chain enzymes were performed 72h post siRNA transfection. Western blots show the efficacy of siRNA-mediated knockdown of the targeted proteins. Image collected and cropped by CiteAb from the following open publication (<https://pubmed.ncbi.nlm.nih.gov/36306280>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Publications

Yasuhara T, Xing Y, Bauer N et al. Condensates induced by transcription inhibition localize active chromatin to nucleoli *Molecular Cell* 2022-06-01 [PMID: 35662392] (Human)

Moghimi S, Viktorova EG, Gabaglio S et al. A Proximity biotinylation assay with a host protein bait reveals multiple factors modulating enterovirus replication *PLOS Pathogens* 2022-10-28 [PMID: 36306280] (Human)

Sikorski K, Mehta A et al. A high-throughput pipeline for validation of antibodies. *Nat Methods* 2018-01-11 [PMID: 30377371] (Human)

Details:

Antibody validation based on denaturing gel electrophoresis of biotinylated cell lysates (PAGE) followed by mass spectrometry (MS) and antibody array analysis (MAP).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-92686

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-92686

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

