

# Product Datasheet

## ARMH1 Recombinant Protein Antigen NBP1-93486PEP

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-93486PEP](http://www.novusbio.com/NBP1-93486PEP)

Updated 11/21/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-93486PEP](http://www.novusbio.com/reviews/destination/NBP1-93486PEP)



**NBP1-93486PEP****ARMH1 Recombinant Protein Antigen**

<b>Product Information</b>	
<b>Unit Size</b>	0.1 ml
<b>Concentration</b>	Please see the vial label for concentration. If unlisted please contact technical services.
<b>Storage</b>	Store at -20C. Avoid freeze-thaw cycles.
<b>Preservative</b>	No Preservative
<b>Purity</b>	>80% by SDS-PAGE and Coomassie blue staining
<b>Buffer</b>	PBS and 1M Urea, pH 7.4.
<b>Target Molecular Weight</b>	25 kDa
<b>Product Description</b>	
<b>Description</b>	<p>A recombinant protein antigen with a N-terminal His6-ABP tag corresponding to human C1orf228.</p> <p><b>Source:</b> <i>E. coli</i></p> <p><b>Amino Acid Sequence:</b> VLGTMHLEVQYEAIELIKDLVGYDVRQALLKGLVALLIPSVKEISKLQAKILSDPS VLQLTPSLPMF</p> <p><b>Fusion Tag:</b> N-terminal His6ABP (ABP = Albumin Binding Protein derived from Streptococcal Protein G)</p> <p><b>This product is intended to be used as a blocking antigen for antibody competition assays. Any other use of this antigen is done at the risk of the user. The use of this product for commercial production is strictly prohibited. Please contact technical support if you have any questions.</b></p>
<b>Gene ID</b>	339541
<b>Gene Symbol</b>	C1orf228
<b>Species</b>	Human
<b>Product Application Details</b>	
<b>Applications</b>	Antibody Competition
<b>Recommended Dilutions</b>	Antibody Competition 10 - 100 molar excess
<b>Application Notes</b>	<p>This recombinant antigen is only intended to be used as a blocking agent to confirm antibody specificity with the corresponding antibody, catalog number NBP1-93486.</p> <p>It is purified by IMAC chromatography, and the expected concentration is greater than 0.5 mg/ml.</p> <p>For current lot information, including availability, please contact our technical support team click <a href="mailto:nb-technical@bio-techne.com">nb-technical@bio-techne.com</a></p>



### Novus Biologicals USA

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### Bio-Techne Canada

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### Bio-Techne Ltd

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### General Contact Information

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### Products Related to NBP1-93486PEP

---

210-TA-005	TNF-alpha [Unconjugated]
NBP1-93486	ARMH1 Antibody - BSA Free
M6000B-1	IL-6 [HRP]
6507-IL-010/CF	IL-4 [Unconjugated]

---

### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-93486PEP](http://www.novusbio.com/reviews/submit/NBP1-93486PEP)

Earn gift cards/discounts by submitting a publication using this product:

[www.novusbio.com/publications](http://www.novusbio.com/publications)

