

Product Datasheet

ARHGAP36 Antibody - BSA Free

NBP1-93894

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-93894

Updated 3/4/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-93894



NBP1-93894

ARHGAP36 Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Host	Rabbit
Gene ID	158763
Gene Symbol	ARHGAP36
Species	Human
Reactivity Notes	Immunogen displays the following percentage of sequence identity for non-tested species: Mouse (81%).
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: NWDVLFQVPPHIQRQVAKRVWKSSPEALDFIRRRNLRKIQSARIKMEEDALLSD PVETSAEARAALVAQSKPSDEGSSEEPVPSGTARSHDDEEGAGNPPIPEQD

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:50 - 1:200
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. Immunocytochemistry/Immunofluorescence Fixation Permeabilization: Use PFA/Triton X-100.

Images


Western Blot: ARHGAP36 Antibody [NBP1-93894] - Analysis in human cell line SH-SY5Y.




Immunohistochemistry-Paraffin: ARHGAP36 Antibody [NBP1-93894] - Staining of human skin shows no positivity in squamous epithelial cells as expected.




Immunohistochemistry-Paraffin: ARHGAP36 Antibody [NBP1-93894] - Staining of human adrenal gland shows high expression.




Immunohistochemistry-Paraffin: ARHGAP36 Antibody [NBP1-93894] - Staining of human pancreas shows low expression as expected.




Immunohistochemistry-Paraffin: ARHGAP36 Antibody [NBP1-93894] - Staining of human pituitary gland shows strong cytoplasmic positivity in cells in anterior lobe.




Immunohistochemistry-Paraffin: ARHGAP36 Antibody [NBP1-93894] - Staining in human adrenal gland and pancreas tissues using anti-ARHGAP36 antibody. Corresponding ARHGAP36 RNA-seq data are presented for the same tissues.




Immunohistochemistry-Paraffin: ARHGAP36 Antibody [NBP1-93894] - Staining of human adrenal gland shows strong cytoplasmic positivity in glandular cells.




Immunohistochemistry-Paraffin: ARHGAP36 Antibody [NBP1-93894] - Staining of human kidney shows no positivity in cells in tubules as expected.




Immunohistochemistry-Paraffin: ARHGAP36 Antibody [NBP1-93894] - Staining of human pancreas shows no positivity in exocrine glandular cells as expected.




Immunohistochemistry-Paraffin: ARHGAP36 Antibody [NBP1-93894] - Staining of human pituitary gland shows strong cytoplasmic positivity in cells in anterior lobe.



Analysis in human cell line SH-SY5Y.



Analysis in human adrenal gland and kidney tissues using NBP1-93894 antibody. Corresponding ARHGAP36 RNA-seq data are presented for the same tissues.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-93894

NBP1-93894PEP	ARHGAP36 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-93894

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

