

Product Datasheet

DYKDDDDK Epitope Tag Antibody (29E4.G7) NBP1-97410-0.1mg

Unit Size: 0.1 mg

Store at -20C short term. Aliquot and store at -80C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 3 Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-97410

Updated 11/11/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-97410



NBP1-97410-0.1mg

DYKDDDDK Epitope Tag Antibody (29E4.G7)

Product Information	
Unit Size	0.1 mg
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C short term. Aliquot and store at -80C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	29E4.G7
Preservative	0.01% Sodium Azide
Isotype	IgG2a Kappa
Purity	Protein A purified
Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Product Description	
Description	<p>This product is an IgG fraction antibody purified from ascites by Protein A chromatography followed by extensive dialysis against the buffer stated above. The purified antibody is directed against the FLAG(TM) motif and is useful in determining its presence in various assays where the epitope tag is present at either the amino or carboxy terminus of recombinant proteins</p> <p>Store vial at -20C prior to opening. Aliquot contents and freeze at -20C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4C as an undiluted liquid. Dilute only prior to immediate use.</p>
Host	Mouse
Species	Epitope Tag
Specificity/Sensitivity	The purified antibody is directed against the FLAG(TM) motif and is useful in determining its presence in various assays where the epitope tag is present at either the amino or carboxy terminus of recombinant proteins. This monoclonal anti-FLAG(TM) tag antibody detects over-expressed proteins containing the FLAG(TM) epitope tag. In western blotting of bacterial extracts, the antibody does not cross-react with endogenous proteins.
Immunogen	This DYKDDDDK Epitope Tag Antibody (29E4.G7) was produced in mice by repeated immunizations with a synthetic peptide corresponding to the FLAG(TM) epitope tag peptide DYKDDDDK Epitope Tag Epitope Tag (Asp-Tyr-Lys-Asp-Asp-Asp-Lys) conjugated to KLH using maleimide. Residues of glycine and cysteine were added to the termini to facilitate coupling.
Notes	<p>FLAG(R) and ANTI-FLAG(R) are registered trademarks of Sigma-Aldrich Co. LLC.</p> <p>FLAG(R) and ANTI-FLAG(R) are registered trademarks of Sigma-Aldrich Co. LLC.</p>

Product Application Details	
Applications	Western Blot, Simple Western, ELISA, Flow Cytometry, Immunohistochemistry
Recommended Dilutions	Western Blot 1:2000-1:10000, Simple Western, Flow Cytometry, ELISA 1:150000-1:250000, Immunohistochemistry 1:1000-1:5000



Application Notes

This product has been tested by ELISA and western blot and is optimally suited for monitoring the expression of FLAG(TM) tagged fusion proteins. As such, this antibody can be used to identify fusion proteins containing the FLAG(TM) epitope. The antibody recognizes the epitope tag fused to either the amino- or carboxy- termini of targeted proteins. The epitope tag peptide sequence was first derived from the 11-amino-acid leader peptide of the gene-10 product from bacteriophage T7. DYKDDDDK is the most commonly used hydrophilic octapeptide tag.


Use in Western Blot and Simple Western reported in a Verified Customer Review. Dilution for use in Flow Cytometry is user optimized.

See [Simple Western Antibody Database](#) for Simple Western validation: separated by Size



Images


Western Blot: DYKDDDDK Epitope Tag Antibody (29E4.G7) [NBP1-97410] - Use of DYKDDDDK Epitope Tag Antibody (29E4.G7) [NBP1-97410] in Western Blot. Flag expression is significantly increased in 293 cells transfected with a flag-fusion protein. Image from a verified customer review.




Western Blot: DYKDDDDK Epitope Tag Antibody (29E4.G7) [NBP1-97410] - Twenty-four clones were randomly selected and grown up from glycerol stocks by inoculating 0.5 mL 2X YT medium. Expression of recombinant proteins was induced by the addition of IPTG. Proteins were purified by nickel affinity chromatography and eluted in 40 uL. Samples were diluted 10-fold, transferred to nitrocellulose membrane and blotted using Mab DYKDDDDK Epitope Tag Antibody (29E4.G7) [NBP1-97410].



Western Blot: DYKDDDDK Epitope Tag Antibody (29E4.G7) [NBP1-97410] - Antibody detects both C-terminal linked and N-terminal linked FLAG tagged recombinant proteins.



ELISA: DYKDDDDK Epitope Tag Antibody (29E4.G7) [NBP1-97410] - Detection of DYKDDDDK conjugated proteins using the DyLight 405 conjugated version of DYKDDDDK Epitope Tag Antibody (29E4.G7).



Emission	Color	DyLight™ Dye	Ex/Em (nm)	ϵ ($M^{-1}cm^{-1}$)	Similar Dyes
Blue		405	400/420	30,000	Alexa™ 405, Cascade Blue
Green		498	493/516	70,000	Alexa™ 498, Cy2®, FITC
Yellow		549	550/568	150,000	Alexa™ 546, Alexa 555, Cy3®, TRITC
Red		649	646/674	250,000	Alexa™ 647, Cy5®
Near Infrared		680	682/715	140,000	Alexa™ 680, Cy5.5®, IRDye™ 700
Infrared		800	770/794	270,000	IRDye™ 800

Publications

Zhang T, Cui S, Xiong X et al. PIH1D3-knockout rats exhibit full ciliopathy features and dysfunctional pre-assembly and loading of dynein arms in motile cilia *Frontiers in cell and developmental biology* 2023-10-12 [PMID: 37900281] (WB)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-97410-0.1mg

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-97410

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

