

Product Datasheet

COASY Antibody (OTI3B3)

NBP2-01081

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-01081

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-01081



NBP2-01081

COASY Antibody (OTI3B3)

Product Information	
Unit Size	0.1 ml
Concentration	0.75 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI3B3
Preservative	0.02% Sodium Azide
Isotype	IgG2a
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.3), 1.0% BSA and 50% Glycerol
Target Molecular Weight	62.1 kDa


Product Description	
Description	Novus Biologicals Mouse COASY Antibody (OTI3B3) (NBP2-01081) is a monoclonal antibody validated for use in IHC, WB, Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	80347
Gene Symbol	COASY
Species	Human, Mouse, Rat
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions.
Immunogen	Full length human recombinant protein of human COASY(NP_079509) produced in HEK293T cell.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:2000, Flow Cytometry 1:100, Immunohistochemistry 1:150, Immunocytochemistry/ Immunofluorescence 1:100, Immunohistochemistry-Paraffin 1:150




Images


Western Blot: COASY Antibody (OTI3B3) [NBP2-01081] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY COASY (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-COASY.




Immunocytochemistry/Immunofluorescence: COASY Antibody (OTI3B3) [NBP2-01081] - Staining of COS7 cells transiently transfected by pCMV6-ENTRY COASY.




Immunohistochemistry-Paraffin: COASY Antibody (OTI3B3) [NBP2-01081] - Staining of paraffin-embedded Human tonsil using anti-COASY mouse monoclonal antibody.




Flow Cytometry: COASY Antibody (OTI3B3) [NBP2-01081] - HEK293T cells transfected with either COASY Human Tagged ORF Clone overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-COASY antibody, and then analyzed by flow cytometry.




Immunohistochemistry-Paraffin: COASY Antibody (OT13B3) [NBP2-01081] - Staining of paraffin-embedded Adenocarcinoma of Human endometrium tissue using anti-COASY mouse monoclonal antibody.



Immunohistochemistry-Paraffin: COASY Antibody (OT13B3) [NBP2-01081] - Staining of paraffin-embedded Adenocarcinoma of Human ovary tissue using anti-COASY mouse monoclonal antibody.




Immunohistochemistry-Paraffin: COASY Antibody (OT13B3) [NBP2-01081] - Staining of paraffin-embedded Carcinoma of Human lung tissue using anti-COASY mouse monoclonal antibody.




Immunohistochemistry-Paraffin: COASY Antibody (OT13B3) [NBP2-01081] - Staining of paraffin-embedded Human breast tissue using anti-COASY mouse monoclonal antibody.



Immunohistochemistry-Paraffin: COASY Antibody (OT13B3) [NBP2-01081] - Staining of paraffin-embedded Human Kidney tissue using anti-COASY mouse monoclonal antibody.



Immunohistochemistry-Paraffin: COASY Antibody (OT13B3) [NBP2-01081] - Staining of paraffin-embedded Human pancreas tissue using anti-COASY mouse monoclonal antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-01081

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-01081

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

