

Product Datasheet

FMRP Antibody (OT11C6) NBP2-01770

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-01770

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-01770



NBP2-01770

FMRP Antibody (OT11C6)

Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OT11C6
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.3), 1.0% BSA and 50% Glycerol
Target Molecular Weight	71 kDa


Product Description	
Description	Novus Biologicals Mouse FMRP Antibody (OT11C6) (NBP2-01770) is a monoclonal antibody validated for use in IHC, WB, Flow and ICC/IF. Anti-FMRP Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	2332
Gene Symbol	FMR1
Species	Human, Mouse, Rat, Canine, Monkey
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions.
Immunogen	Human recombinant protein fragment corresponding to amino acids 36-279 of human FMR1(NP_002015) produced in E.coli.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-2000, Flow Cytometry 1:100, Immunohistochemistry 1:150, Immunocytochemistry/ Immunofluorescence 1:100, Immunohistochemistry-Paraffin 1:150




Images


Western Blot: FMRP Antibody (1C6) [NBP2-01770] Analysis of extracts (35ug) from 9 different cell lines by using anti-FMRP monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).




Immunocytochemistry/Immunofluorescence: FMRP Antibody (1C6) [NBP2-01770] - Staining of COS7 cells transiently transfected by pCMV6-ENTRY FMRP.




Immunohistochemistry-Paraffin: FMRP Antibody (1C6) [NBP2-01770] - Staining of paraffin-embedded Human pancreas tissue using anti-FMRP mouse monoclonal antibody.




Flow Cytometry: FMRP Antibody (1C6) [NBP2-01770] - Analysis of Jurkat cells, using anti-FMRP antibody, (Red), compared to a nonspecific negative control antibody (Blue).




Western Blot: FMRP Antibody (1C6) [NBP2-01770] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY FMRP (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FMRP.




Immunohistochemistry-Paraffin: FMRP Antibody (1C6) [NBP2-01770] - Staining of paraffin-embedded Adenocarcinoma of Human breast tissue using anti-FMRP mouse monoclonal antibody.




Immunohistochemistry-Paraffin: FMRP Antibody (1C6) [NBP2-01770] - Staining of paraffin-embedded Adenocarcinoma of Human colon tissue using anti-FMRP mouse monoclonal antibody.



Immunohistochemistry-Paraffin: FMRP Antibody (1C6) [NBP2-01770] - Staining of paraffin-embedded Carcinoma of Human kidney tissue using anti-FMRP mouse monoclonal antibody.



Immunohistochemistry-Paraffin: FMRP Antibody (1C6) [NBP2-01770] - Staining of paraffin-embedded Human colon tissue using anti-FMRP mouse monoclonal antibody.



Publications

Yildirim Z, Baboo S, Hamid SM et al. Intercepting IRE1 kinase-FMRP signaling prevents atherosclerosis progression EMBO molecular medicine 2022-02-22 [PMID: 35191199] (WB)

Yang YM, Arsenault J, Bah A et al. Identification of a molecular locus for normalizing dysregulated GABA release from interneurons in the Fragile X brain. Mol. Psychiatry. 2018-09-17 [PMID: 30224722] (WB, Mouse)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-01770

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-01770

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

