

# Product Datasheet

## PSMB7 Antibody (OT11E1)

### NBP2-01832

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

#### Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-01832](http://www.novusbio.com/NBP2-01832)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-01832](http://www.novusbio.com/reviews/destination/NBP2-01832)



**NBP2-01832**

PSMB7 Antibody (OT11E1)


Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OT11E1
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.3), 1.0% BSA and 50% Glycerol
Target Molecular Weight	29.8 kDa

Product Description	
Description	Novus Biologicals Mouse PSMB7 Antibody (OT11E1) (NBP2-01832) is a monoclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	5695
Gene Symbol	PSMB7
Species	Human, Mouse, Rat, Monkey, Primate
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions.
Immunogen	Human recombinant protein fragment corresponding to amino acids 58-277 of human PSMB7(NP_002790) produced in E.coli.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-2000, Immunohistochemistry 1:150, Immunohistochemistry-Paraffin 1:150

**Images**


Western Blot: PSMB7 Antibody (OT11E1) [NBP2-01832] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY Proteasome 20S beta 7 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-Proteasome 20S beta 7.




Immunohistochemistry-Paraffin: PSMB7 Antibody (OT11E1) [NBP2-01832] - Staining of paraffin-embedded Human prostate tissue using anti-Proteasome 20S beta 7 mouse monoclonal antibody.




Western Blot: PSMB7 Antibody (OT11E1) [NBP2-01832] - Analysis of extracts (35ug) from 9 different cell lines by using anti-Proteasome 20S beta 7 monoclonal antibody.




Immunohistochemistry-Paraffin: PSMB7 Antibody (OT11E1) [NBP2-01832] - Staining of paraffin-embedded Adenocarcinoma of Human breast tissue using anti-Proteasome 20S beta 7 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: PSMB7 Antibody (OT11E1) [NBP2-01832] - Staining of paraffin-embedded Adenocarcinoma of Human colon tissue using anti-Proteasome 20S beta 7 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: PSMB7 Antibody (OT11E1) [NBP2-01832] - Staining of paraffin-embedded Carcinoma of Human bladder tissue using anti-Proteasome 20S beta 7 mouse monoclonal antibody.



Immunohistochemistry-Paraffin: PSMB7 Antibody (OT11E1) [NBP2-01832] - Staining of paraffin-embedded Human breast tissue using anti-Proteasome 20S beta 7 mouse monoclonal antibody.



Immunohistochemistry-Paraffin: PSMB7 Antibody (OT11E1) [NBP2-01832] - Staining of paraffin-embedded Human colon tissue using anti-Proteasome 20S beta 7 mouse monoclonal antibody.



## Publications

Halbert D, Domenyuk V, Spetzler D et al. Aptamers and uses thereof United States Patent Application US 9958448 B2 2018-01-01



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-01832**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-01832](http://www.novusbio.com/reviews/submit/NBP2-01832)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

