

# Product Datasheet

## Ketohexokinase Antibody (OTI1D8) NBP2-02639

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-02639](http://www.novusbio.com/NBP2-02639)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-02639](http://www.novusbio.com/reviews/destination/NBP2-02639)



**NBP2-02639****Ketohexokinase Antibody (OT11D8)**


<b>Product Information</b>	
<b>Unit Size</b>	0.1 ml
<b>Concentration</b>	0.57 mg/ml
<b>Storage</b>	Store at -20C. Avoid freeze-thaw cycles.
<b>Clonality</b>	Monoclonal
<b>Clone</b>	OT11D8
<b>Preservative</b>	0.02% Sodium Azide
<b>Isotype</b>	IgG1
<b>Purity</b>	Immunogen affinity purified
<b>Buffer</b>	PBS (pH 7.3), 1.0% BSA and 50% Glycerol
<b>Target Molecular Weight</b>	32.5 kDa

<b>Product Description</b>	
<b>Description</b>	Novus Biologicals Mouse Ketohexokinase Antibody (OT11D8) (NBP2-02639) is a monoclonal antibody validated for use in IHC, WB, Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
<b>Host</b>	Mouse
<b>Gene ID</b>	3795
<b>Gene Symbol</b>	KHK
<b>Species</b>	Human
<b>Immunogen</b>	Full length human recombinant protein of human KHK (NP_000212) produced in HEK293T cell.


<b>Product Application Details</b>	
<b>Applications</b>	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
<b>Recommended Dilutions</b>	Western Blot 1:2000, Flow Cytometry 1:100, Immunohistochemistry 1:50, Immunocytochemistry/ Immunofluorescence 1:100, Immunohistochemistry-Paraffin 1:50

**Images**


Western Blot: ketohexokinase Antibody (1D8) [NBP2-02639] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ketohexokinase (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ketohexokinase.




Immunocytochemistry/Immunofluorescence: Ketohekinase Antibody (OT11D8) [NBP2-02639] - Staining of COS7 cells transiently transfected by pCMV6-ENTRY ketohekinase.




Immunohistochemistry-Paraffin: Ketohekinase Antibody (1D8) [NBP2-02639] - ketohekinase Antibody (1D8) [NBP2-02639] - Staining of paraffin-embedded Human thyroid tissue using anti-ketohekinase mouse monoclonal antibody.




Flow Cytometry: Ketohekinase Antibody (1D8) [NBP2-02639] - ketohekinase Antibody (1D8) [NBP2-02639] - HEK293T cells transfected with either overexpression plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-ketohekinase antibody, and then analyzed by flow cytometry.




Immunohistochemistry-Paraffin: Ketohekinase Antibody (1D8) [NBP2-02639] - ketohekinase Antibody (1D8) [NBP2-02639] - Staining of paraffin-embedded Adenocarcinoma of Human colon tissue using anti-ketohekinase mouse monoclonal antibody.




Immunohistochemistry-Paraffin: Ketohekinase Antibody (1D8) [NBP2-02639] - ketohekinase Antibody (1D8) [NBP2-02639] - Staining of paraffin-embedded Adenocarcinoma of Human ovary tissue using anti-ketohekinase mouse monoclonal antibody.



Immunohistochemistry-Paraffin: Ketohekinase Antibody (1D8) [NBP2-02639] - ketohekinase Antibody (1D8) [NBP2-02639] - Staining of paraffin-embedded Human Kidney tissue using anti-ketohekinase mouse monoclonal antibody.



Immunohistochemistry-Paraffin: Ketohekinase Antibody (1D8) [NBP2-02639] - ketohekinase Antibody (1D8) [NBP2-02639] - Staining of paraffin-embedded Human liver tissue using anti-ketohekinase mouse monoclonal antibody.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-02639**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-02639](http://www.novusbio.com/reviews/submit/NBP2-02639)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

