

Product Datasheet

USP36 Antibody (OTI7G3)

NBP2-03830

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-03830

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-03830



NBP2-03830

USP36 Antibody (OTI7G3)


Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI7G3
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.3), 1.0% BSA and 50% Glycerol
Target Molecular Weight	122.7 kDa

Product Description	
Description	Novus Biologicals Mouse USP36 Antibody (OTI7G3) (NBP2-03830) is a monoclonal antibody validated for use in IHC, WB and ICC/IF. Anti-USP36 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	57602
Gene Symbol	USP36
Species	Human, Canine
Marker	Nucleolus Marker
Immunogen	Human recombinant protein fragment corresponding to amino acids 589-972 of human USP36 (NP_079366) produced in E.coli.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Knockdown Validated
Recommended Dilutions	Western Blot 1:500-2000, Immunohistochemistry 1:150, Immunocytochemistry/Immunofluorescence 1:100, Immunohistochemistry-Paraffin 1:150, Knockdown Validated

Images


Western Blot: USP36 Antibody (OTI7G3) [NBP2-03830] - Analysis of extracts (35ug) from 9 different cell lines by using anti-USP36 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).




Immunocytochemistry/Immunofluorescence: USP36 Antibody (OTI7G3) [NBP2-03830] - Staining of COS7 cells transiently transfected by pCMV6-ENTRY USP36.




Immunohistochemistry-Paraffin: USP36 Antibody (OTI7G3) [NBP2-03830] - Staining of paraffin-embedded Human liver tissue using anti-USP36 mouse monoclonal antibody.




Western Blot: USP36 Antibody (OTI7G3) [NBP2-03830] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY USP36 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-USP36.




Immunohistochemistry-Paraffin: USP36 Antibody (OTI7G3) [NBP2-03830] - Staining of paraffin-embedded Carcinoma of Human kidney tissue using anti-USP36 mouse monoclonal antibody.



Immunohistochemistry-Paraffin: USP36 Antibody (OTI7G3) [NBP2-03830] - Staining of paraffin-embedded Carcinoma of Human liver tissue using anti-USP36 mouse monoclonal antibody.



Immunohistochemistry-Paraffin: USP36 Antibody (OTI7G3) [NBP2-03830] - Staining of paraffin-embedded Carcinoma of Human lung tissue using anti-USP36 mouse monoclonal antibody.



Publications

Iden M, Tsaih SW, Huang YW et al. Multi-omics mapping of human papillomavirus integration sites illuminates novel cervical cancer target genes Br J Cancer 2021-09-15 [PMID: 34526665] (WB, KD, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-03830

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-03830

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

