

Product Datasheet

Proteasome 20S alpha 6 Antibody (OTI4C9) NBP2-03836

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-03836

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-03836



NBP2-03836

Proteasome 20S alpha 6 Antibody (OTI4C9)


Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI4C9
Preservative	0.02% Sodium Azide
Isotype	IgG2b
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.3), 1.0% BSA and 50% Glycerol
Target Molecular Weight	27.2 kDa

Product Description	
Description	Novus Biologicals Mouse Proteasome 20S alpha 6 Antibody (OTI4C9) (NBP2-03836) is a monoclonal antibody validated for use in IHC, WB and Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	5687
Gene Symbol	PSMA6
Species	Human, Mouse, Rat, Canine, Primate, Monkey
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions.
Immunogen	Full length human recombinant protein of human PSMA6 (NP_002782) produced in E.coli.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-2000, Flow Cytometry 1:100, Immunohistochemistry 1:150, Immunohistochemistry-Paraffin 1:150

Images


Western Blot: Proteasome 20S alpha 6 Antibody (4C9) [NBP2-03836] Analysis of extracts (35ug) from 9 different cell lines by using anti-Proteasome 20S alpha 6 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).




Immunohistochemistry-Paraffin: Proteasome 20S alpha 6 Antibody (4C9) [NBP2-03836] - Staining of paraffin-embedded Human pancreas tissue using anti-Proteasome 20S alpha 6 mouse monoclonal antibody.




Flow Cytometry: Proteasome 20S alpha 6 Antibody (4C9) [NBP2-03836] - HEK293T cells transfected with either overexpression plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-Proteasome 20S alpha 6 antibody, and then analyzed by flow cytometry.




Western Blot: Proteasome 20S alpha 6 Antibody (4C9) [NBP2-03836] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY Proteasome 20S alpha 6 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-Proteasome 20S alpha 6.




Immunohistochemistry-Paraffin: Proteasome 20S alpha 6 Antibody (4C9) [NBP2-03836] - Staining of paraffin-embedded Adenocarcinoma of Human endometrium tissue using anti-Proteasome 20S alpha 6 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: Proteasome 20S alpha 6 Antibody (4C9) [NBP2-03836] - Staining of paraffin-embedded Adenocarcinoma of Human ovary tissue using anti-Proteasome 20S alpha 6 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: Proteasome 20S alpha 6 Antibody (4C9) [NBP2-03836] - Staining of paraffin-embedded Carcinoma of Human bladder tissue using anti-Proteasome 20S alpha 6 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: Proteasome 20S alpha 6 Antibody (4C9) [NBP2-03836] - Staining of paraffin-embedded Carcinoma of Human kidney tissue using anti-Proteasome 20S alpha 6 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: Proteasome 20S alpha 6 Antibody (4C9) [NBP2-03836] - Staining of paraffin-embedded Carcinoma of Human lung tissue using anti-Proteasome 20S alpha 6 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: Proteasome 20S alpha 6 Antibody (4C9) [NBP2-03836] - Staining of paraffin-embedded Carcinoma of Human prostate tissue using anti-Proteasome 20S alpha 6 mouse monoclonal antibody.



Immunohistochemistry-Paraffin: Proteasome 20S alpha 6 Antibody (4C9) [NBP2-03836] - Staining of paraffin-embedded Carcinoma of Human thyroid tissue using anti-Proteasome 20S alpha 6 mouse monoclonal antibody.



Immunohistochemistry-Paraffin: Proteasome 20S alpha 6 Antibody (4C9) [NBP2-03836] - Staining of paraffin-embedded Human Kidney tissue using anti-Proteasome 20S alpha 6 mouse monoclonal antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-03836

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-03836

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

