

Product Datasheet

NLRP3/NALP3 Antibody (25N10E9) - BSA Free NBP2-03948

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 3

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-03948

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-03948



NBP2-03948

NLRP3/NALP3 Antibody (25N10E9) - BSA Free

| Product Information | |
|-----------------------------|--|
| Unit Size | 0.1 mg |
| Concentration | 1.0 mg/ml |
| Storage | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality | Monoclonal |
| Clone | 25N10E9 |
| Preservative | 0.02% Sodium Azide |
| Isotype | IgG2b Kappa |
| Purity | Protein G purified |
| Buffer | PBS |
| Product Description | |
| Description | Novus Biologicals Mouse NLRP3/NALP3 Antibody (25N10E9) - BSA Free (NBP2-03948) is a monoclonal antibody validated for use in IHC, WB and Simple Western. Anti-NLRP3/NALP3 Antibody: Cited in 3 publications. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Mouse |
| Gene ID | 114548 |
| Gene Symbol | NLRP3 |
| Species | Mouse |
| Specificity/Sensitivity | The epitope that 25N10E9 binds to the NOD region of NLRP3/NALP3 has not been mapped. |
| Immunogen | A recombinant protein corresponding to the NOD domain of murine NLRP3/NALP3 Antibody (25N10E9) (amino acids 216-385) was used as the immunogen for this antibody. |
| Product Application Details | |
| Applications | Western Blot, Simple Western |
| Recommended Dilutions | Western Blot 2-5 ug/ml, Simple Western 1:100 |
| Application Notes | In Simple Western only 10 - 15 uL of the recommended dilution is used per data point. See Simple Western Antibody Database for Simple Western validation: Tested in Bone marrow, separated by Size, antibody dilution of 1:100 |




Images


Simple Western: NLRP3/NALP3 Antibody (25N10E9) [NBP2-03948] - Simple Western lane view shows a specific band for NALP3 in 0.5 mg/ml of RAW lysate. This experiment was performed under reducing conditions using the 12-230 kDa separation system.



Western Blot: NLRP3/NALP3 Antibody (25N10E9) [NBP2-03948] - Analysis of Raw 264 cell lysate using NLRP3/NALP3 antibody (NBP2-03948) at 2 ug/ml.



Western Blot: NLRP3/NALP3 Antibody (25N10E9) [NBP2-03948] - Analysis using the NALP3 monoclonal antibody in mouse NALP3 transfected cell lysate (A), human NALP3 transfected cell lysate (B), mouse RAW lysate (C) and the recombinant mouse NLRP3 NOD domain protein immunogen (D) probed with 2 ug/ml antibody. Goat anti-mouse IgG HRP secondary antibody and PicoTect ECL substrate solution were used for this test.



Publications

Wang X, Cheng W, Chen X et al. Inhibition of CEBPB Attenuates Lupus Nephritis via Regulating Pim-1 Signaling Mediators of Inflammation 2022-09-20 [PMID: 36248187] (Immunohistochemistry, Mouse)

Chao YY, Puhach A, Frieser D et al. Human TH17 cells engage gasdermin E pores to release IL-1 β on NLRP3 inflammasome activation Nature immunology 2023-02-01 [PMID: 36604548] (Simple Western, Human)

Perez-Martinez PI, Rojas-Espinosa O, Hernandez-Chavez VG et al. Anti-inflammatory effect of omega unsaturated fatty acids and dialysable leucocyte extracts on collagen-induced arthritis in DBA/1 mice Int J Exp Pathol 2020-05-28 [PMID: 32459025] (IF/IHC, Mouse)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-03948

| | |
|------------------|--|
| NBP2-33376H | Blue Marker Antibody (6F4-F6) [HRP] |
| HAF007 | Goat anti-Mouse IgG Secondary Antibody [HRP] |
| NB7539 | Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP] |
| NBP1-43317-0.5mg | Mouse IgG2b Kappa Light Chain Isotype Control (MG2b) |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-03948

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

