

Product Datasheet

ILT5/CD85a/LILRB3 Antibody (MM0413-9S32) - Azide and BSA Free NBP2-11729

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-11729

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-11729



NBP2-11729

ILT5/CD85a/LILRB3 Antibody (MM0413-9S32) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	LYOPH mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	MM0413-9S32
Preservative	No Preservative
Reconstitution Instructions	Reconstitute with sterilized PBS to a final concentration of 0.5 mg/ml.
Isotype	IgG2a
Purity	Protein G purified
Buffer	Lyophilized from a 0.2 um filtered solution in PBS. 0.025 mg size is provided in liquid form, PBS


Product Description	
Description	Novus Biologicals Mouse ILT5/CD85a/LILRB3 Antibody (MM0413-9S32) - Azide and BSA Free (NBP2-11729) is a monoclonal antibody validated for use in IHC, WB and Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	11025
Gene Symbol	LILRB3
Species	Human
Immunogen	Human recombinant ILT5/CD85a/LILRB3

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunohistochemistry, Block/Neutralize, CyTOF-ready
Recommended Dilutions	Western Blot 1:500-1:2000, Flow Cytometry 1:10-1:1000, Immunohistochemistry 1:50-1:100, Immunohistochemistry-Paraffin, CyTOF-ready, Block/Neutralize
Application Notes	This antibody is Cytof ready.




Images

Immunohistochemistry-Paraffin: ILT5/CD85a/LILRB3 Antibody (MM0413-9S32) [NBP2-11729] - 10% Buffer formalin fixed and paraffine embedded human placental tissue section (4um) is subjected IHC staining using NBP2-11729. Tissue section was pretreated in citric buffer (pH 6.0) with microwave for antigen retrieval before IHC is applied.



Flow Cytometry: ILT5/CD85a/LILRB3 Antibody (MM0413-9S32) [NBP2-11729] - Analysis using Allophycocyanin conjugate of NBP2-11729. A cell surface stain was performed on hPBMCs with ILT5 antibody (MM0413-9S32) NBP2-11729 (top image) and a matched isotype control NBP2-27287 (bottom image). Cells were incubated in an antibody dilution of 1:200 for 20 minutes at room temperature. Both antibodies were conjugated to APC. A co-stain was performed using CD14 antibody (M5E2) NB100-77758PE.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-11729

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-11729

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

