

Product Datasheet

TGF-beta RIII Antibody (MM0057-5G9) - Azide and BSA Free NBP2-11920

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-11920

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-11920



NBP2-11920

TGF-beta RIII Antibody (MM0057-5G9) - Azide and BSA Free


| Product Information | |
|-----------------------------|---|
| Unit Size | 0.1 mg |
| Concentration | LYOPH mg/ml |
| Storage | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality | Monoclonal |
| Clone | MM0057-5G9 |
| Preservative | No Preservative |
| Reconstitution Instructions | Reconstitute with sterilized PBS to a final concentration of 0.2 mg/ml. |
| Isotype | IgG1 |
| Purity | Protein G purified |
| Buffer | Lyophilized from a 0.2 um filtered solution in PBS. 0.025 mg size is provided in liquid form, PBS |

| Product Description | |
|-------------------------|---|
| Description | Novus Biologicals Mouse TGF-beta RIII Antibody (MM0057-5G9) - Azide and BSA Free (NBP2-11920) is a monoclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Mouse |
| Gene ID | 7049 |
| Gene Symbol | TGFBR3 |
| Species | Human |
| Specificity/Sensitivity | In direct ELISAs, this antibody does not cross-react with recombinant human TGF-beta RI, and TGF-beta RII. |
| Immunogen | Human recombinant TGF-beta RIII (extracellular domain) |

| Product Application Details | |
|-----------------------------|--|
| Applications | Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry |
| Recommended Dilutions | Western Blot 1:100-1:1000, Immunohistochemistry 1:10-1:500, Immunohistochemistry-Paraffin 1:50-1:100 |

Images

Immunohistochemistry-Paraffin: TGF beta Receptor III Antibody (MM0057-5G9) [NBP2-11920] - Staining of Human placenta.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-11920

| | |
|------------------|--|
| NBP2-33376H | Blue Marker Antibody (6F4-F6) [HRP] |
| HAF007 | Goat anti-Mouse IgG Secondary Antibody [HRP] |
| NB7539 | Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP] |
| NBP1-97005-0.5mg | Mouse IgG1 Isotype Control (MG1) |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-11920

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

