

Product Datasheet

IDH3B Antibody - BSA Free

NBP2-14114

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-14114

Updated 12/2/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-14114



NBP2-14114

IDH3B Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Description	Novus Biologicals Rabbit IDH3B Antibody - BSA Free (NBP2-14114) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-IDH3B Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3420
Gene Symbol	IDH3B
Species	Human
Immunogen	This antibody was developed against a recombinant protein corresponding to the amino acids: VKEVFKAAAVPVEFQEHLSEVQNMASEEKLEQVLSSMKENKVAIIGKIHTPME YKGEASYDMRLRRKLDLDFANVVHVKSLPGYMTRHNN


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:2500 - 1:5000, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:2500-1:5000
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100.

Images


Western Blot: IDH3B Antibody [NBP2-14114] - Analysis using Anti-IDH3B antibody NBP2-14114 (A) shows similar pattern to independent antibody NBP2-58706 (B).




Immunocytochemistry/Immunofluorescence: IDH3B Antibody [NBP2-14114] - Staining of human cell line A549 shows localization to mitochondria. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: IDH3B Antibody [NBP2-14114] - Staining of human testis shows strong granular cytoplasmic positivity in cells in seminiferous ducts.




Immunohistochemistry-Paraffin: IDH3B Antibody [NBP2-14114] - Staining of human kidney shows strong granular cytoplasmic positivity in cells in tubules.




Immunohistochemistry-Paraffin: IDH3B Antibody [NBP2-14114] - Staining of human pancreas shows strong granular cytoplasmic positivity in exocrine glandular cells.



Immunohistochemistry-Paraffin: IDH3B Antibody [NBP2-14114] - Staining of human parathyroid gland shows strong granular cytoplasmic positivity in glandular cells.



Immunohistochemistry-Paraffin: IDH3B Antibody [NBP2-14114] - Staining of human placenta shows strong granular cytoplasmic positivity in trophoblastic cells.



Publications

Ruan L, McNamara JT, Zhang X et al. Solid-phase inclusion as a mechanism for regulating unfolded proteins in the mitochondrial matrix *Sci Adv* 2020-08-01 [PMID: 32821848]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-14114

NBP2-14114PEP	IDH3B Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-14114

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

