

Product Datasheet

ACADVL Antibody NBP2-15238

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-15238

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-15238



NBP2-15238

ACADVL Antibody


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS, 1% BSA, 20% Glycerol
Target Molecular Weight	70 kDa

Product Description	
Description	Novus Biologicals Knockout (KO) Validated Rabbit ACADVL Antibody (NBP2-15238) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-ACADVL Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	37
Gene Symbol	ACADVL
Species	Human, Mouse, Rat
Reactivity Notes	Immunogen displays the following percentage of sequence identity for non-tested species: Zebrafish (81%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human ACADVL. The exact sequence is proprietary.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Knockout Validated
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000, Knockout Validated

Images


Immunocytochemistry/Immunofluorescence: ACADVL Antibody [NBP2-15238] - ACADVL antibody detects ACADVL protein at mitochondria by immunofluorescent analysis. Sample: HeLa cells were fixed in ice-cold MeOH for 5 min. Green: ACADVL stained by ACADVL antibody diluted at 1:500. Blue: Fluoroshield with DAPI. Scale bar= 10 um.




Immunohistochemistry-Paraffin: ACADVL Antibody [NBP2-15238] - Rat liver. ACADVL stained by ACADVL antibody [N1C1] diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min.




Western Blot: ACADVL Antibody [NBP2-15238] - Sample (50 ug of whole cell lysate) A: mouse liver 7. 5% SDS PAGE, antibody diluted at 1:3000.




Immunohistochemistry-Paraffin: ACADVL Antibody [NBP2-15238] - H358 xenograft, using ACADVL antibody at 1:250 dilution. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min.




Immunohistochemistry-Paraffin: ACADVL Antibody [NBP2-15238] - Human A549 xenograft. ACADVL antibody [N1C1] diluted at 1:500. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min




Immunohistochemistry-Paraffin: ACADVL Antibody [NBP2-15238] - Mouse kidney. ACADVL stained by ACADVL antibody [N1C1] diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min.




Various whole cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with ACADVL antibody [N1C1] (NBP2-15238) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody. Corresponding RNA expression data for the same cell lines are based on Human Protein Atlas program.



Wild-type (WT) and ACADVL knockout (KO) 293T cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with ACADVL antibody [N1C1] (NBP2-15238) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.



Various whole cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with ACADVL antibody [N1C1] (NBP2-15238) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.



Publications

Marumo T, Maduka CV, Ural E et al. Flavinated SDHA underlies the change in intrinsic optical properties of oral cancers Communications biology 2023-11-09 [PMID: 37945749] (ICC/IF, WB, Human)

Details:

1:500 dilution WB

Tcheng M, Roma A, Ahmed N et al. Very long chain fatty acid metabolism is required in acute myeloid leukemia Blood 2021-03-15 [PMID: 33720355]

Details:

Antibody used in co-immunoprecipitation





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-15238

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-15238

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

