

# Product Datasheet

## Histone H3 [Trimethyl Lys27] Antibody NBP2-16840

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

**Reviews: 1 Publications: 2**

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-16840](http://www.novusbio.com/NBP2-16840)

Updated 9/25/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-16840](http://www.novusbio.com/reviews/destination/NBP2-16840)




**NBP2-16840**

## Histone H3 [Trimethyl Lys27] Antibody


<b>Product Information</b>	
<b>Unit Size</b>	0.1 ml
<b>Concentration</b>	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
<b>Storage</b>	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
<b>Clonality</b>	Polyclonal
<b>Preservative</b>	0.025% Proclin 300
<b>Isotype</b>	IgG
<b>Purity</b>	Antigen Affinity-purified
<b>Buffer</b>	PBS, 1% BSA, 20% Glycerol
<b>Target Molecular Weight</b>	15 kDa
<b>Product Description</b>	
<b>Description</b>	Novus Biologicals Rabbit Histone H3 [Trimethyl Lys27] Antibody (NBP2-16840) is a polyclonal antibody validated for use in IHC, WB, ICC/IF and ChIP. Anti-Histone H3 Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
<b>Host</b>	Rabbit
<b>Gene ID</b>	126961
<b>Gene Symbol</b>	H3C14
<b>Species</b>	Human, Mouse, Rat, Zebrafish
<b>Immunogen</b>	Carrier-protein conjugated synthetic peptide surrounding tri-methyl Lys27 of human Histone H3. The exact sequence is proprietary.
<b>Product Application Details</b>	
<b>Applications</b>	Western Blot, Immunohistochemistry-Paraffin, Dot Blot, Electron Microscopy, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Chromatin Immunoprecipitation (ChIP), Immunohistochemistry Whole-Mount
<b>Recommended Dilutions</b>	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000, Immunohistochemistry-Frozen 1:100-1:1000, Dot Blot 1:100 - 1:2000, Electron Microscopy 1:10 - 1:500, Immunohistochemistry Whole-Mount 1:100-1:500, Chromatin Immunoprecipitation (ChIP) Reported in scientific literature (PMID: 30143902)

## Images


Immunocytochemistry/Immunofluorescence: Histone H3 [Trimethyl Lys27] Antibody [NBP2-16840] - 293T cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: Histone H3K27me3 (Tri-methyl Lys27) stained by Histone H3K27me3 (Tri-methyl Lys27) antibody diluted at 1:500. Red: alpha Tubulin, a cytoskeleton marker, stained by alpha Tubulin antibody [114] diluted at 1:1000.




Immunohistochemistry-Paraffin: Histone H3 [Trimethyl Lys27] Antibody [NBP2-16840] - Mouse intestine. Histone H3K27me3 (Tri-methyl Lys27) stained by Histone H3K27me3 (Tri-methyl Lys27) antibody diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min.




Western Blot: Histone H3 [Trimethyl Lys27] Antibody [NBP2-16840] - Histone H3 [Trimethyl Lys27] Antibody U87 cells were grown under normoxic (DMEM 10% FBS, 16% O<sub>2</sub>), hypoxic (DMEM 10% FBS, 1% O<sub>2</sub>), and neurosphere conditions (DMEM/F12, B-27 supplement, growth factor (10ng/ml FGF and 20ng/ml EGF)). Cell lysate were Western blotted for H3K27me3 and Histone-H3 (loading control).




Immunocytochemistry/Immunofluorescence: Histone H3 [Trimethyl Lys27] Antibody [NBP2-16840] - Histone H3 [Trimethyl Lys27] Antibody Confocal immunofluorescence analysis of paraformaldehyde-fixed A431, using HIST1H3A (tri-Methyl K27) antibody (Green) at 1:500 dilution. Alpha-tubulin filaments are labeled with Alpha-tubulin antibody (Red) at 1:2000.




Immunohistochemistry-Paraffin: Histone H3 [Trimethyl Lys27] Antibody [NBP2-16840] - Hela xenograft, using Histon H3 (tri-Methyl K27) antibody at 1:500 dilution. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min.




Western Blot: Histone H3 [Trimethyl Lys27] Antibody [NBP2-16840] - Various whole cell extracts (30 ug) were separated by 15% SDS-PAGE, and the membrane was blotted with Histone H3 antibody (NBP2-16840) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.




Western Blot: Histone H3 [Trimethyl Lys27] Antibody [NBP2-16840] - Various whole cell extracts (30 ug) were separated by 15% SDS-PAGE, and the membrane was blotted with Histone H3 antibody (NBP2-16840) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




Immunocytochemistry/ Immunofluorescence: Histone H3 [Trimethyl Lys27] Antibody [NBP2-16840] - Histone H3 antibody detects Histone H3 protein at nucleus by immunofluorescent analysis.  
Sample: 293T cells were fixed in 4% paraformaldehyde at RT for 15 min.  
Green: Histone H3 stained by Histone H3 antibody (NBP2-16840) diluted at 1:500.  
Red: alpha Tubulin, a cytoskeleton marker, stained by alpha Tubulin antibody [GT114] diluted at 1:1000.




Western Blot: Histone H3 [Trimethyl Lys27] Antibody [NBP2-16840] - Various whole cell extracts (30 ug) were separated by 15% SDS-PAGE, and the membrane was blotted with Histone H3 antibody (NBP2-16840) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.



Various peptide samples were spotted onto positively charged nylon membrane, and blotted with Histone H3 (tri-Methyl Lys27) antibody (NBP2-16840) diluted at 1:2500.



Various whole cell extracts (30 ug) were separated by 15% SDS-PAGE, and the membrane was blotted with Histone H3K27me3 (trimethyl Lys27) antibody (NBP2-16840) diluted at 1:1000.




Immunofluorescence photomicrographs of frozen sections of mouse liver.


Red: Histone H3K27me3 (trimethyl Lys27) antibody (NBP2-16840) diluted at 1:200. The signal was developed using goat anti-rabbit IgG antibody (Dylight594) .

Blue: Nuclear staining with Hoechst 33342.

br>Antigen Retrieval: Trilogy(TM) (EDTA based, pH 8.0) buffer, 15min



Histone H3K27me3 (trimethyl Lys27) antibody detects Histone H3K27me3 (trimethyl Lys27) protein on zebrafish by whole mount immunohistochemical analysis. Sample: Paraformaldehyde-fixed 2 day-post-fertilization zebrafish embryo. Histone H3K27me3 (trimethyl Lys27) antibody (NBP2-16840) dilution: 1:100.



## Publications

Shen Y, Zhao H, Zhang L Et Al. The roles of DNA methylation and hydroxymethylation at short interspersed nuclear elements in the hypothalamic arcuate nucleus during puberty Molecular Therapy - Nucleic Acids 2021-07-01 [PMID: 34513307] (Chemotaxis)

Aroni MAT, de Oliveira GJPL, Spolidorio LC et al. Loading deproteinized bovine bone with strontium enhances bone regeneration in rat calvarial critical size defects Clin Oral Investig 2018-08-25 [PMID: 30143902] (Chemotaxis, Mouse)



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-16840**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-16840](http://www.novusbio.com/reviews/submit/NBP2-16840)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

