

Product Datasheet

SGK1 Antibody NBP2-19010

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-19010

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-19010



NBP2-19010

SGK1 Antibody


Product Information	
Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	PBS with 0.05% BSA

Product Description	
Description	Novus Biologicals Rabbit SGK1 Antibody (NBP2-19010) is a polyclonal antibody validated for use in IHC and WB. Anti-SGK1 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	6446
Gene Symbol	SGK1
Species	Human
Reactivity Notes	Immunogen displays the following percentage of sequence identity for non-tested species: 100% homologous in new world monkey, macaque, rhesus monkey, opossum, gibbon, orangutan, and chimpanzee; 95% in panda, carolina anole, cow, dog, chinese hamster, zebrafish, horse, mole rat, elephant, platypus, and pig; 90% in mouse, rat, zebra finch, frog, cichlid fish, rabbit, guinea pig, chicken, and wild turkey; 85% in salmon.
Immunogen	A portion of amino acids 360-410 from human SGK1 was used as the immunogen for this antibody. NP_005618.2.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1-3 ug/ml, Immunohistochemistry 1:10-1:500, Immunohistochemistry-Paraffin 5 ug/ml

Images

Western Blot: SGK1 Antibody [NBP2-19010] - Analysis of SGK1 in human lung lysate in the 1) absence and 2) presence of immunizing peptide, 3) human kidney lysate probed with 1 ug/ml of SGK1 antibody. Goat anti-rabbit Ig HRP secondary antibody and PicoTect ECL substrate solution were used for this western blot.



Immunohistochemistry-Paraffin: SGK1 Antibody [NBP2-19010] - analysis of SGK1 in human lung probed with SGK1 antibody at 5 ug/ml.



Publications

Wen W, Wan Z, Ren K et al. Potassium supplementation inhibits IL-17A production induced by salt loading in human T lymphocytes via p38/MAPK-SGK1 pathway. *Experimental and Molecular Pathology*. 2016-03-25 [PMID: 27020669] (WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-19010

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-19010

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

