

# Product Datasheet

## OCT1 Antibody NBP2-21584

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

### Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-21584](http://www.novusbio.com/NBP2-21584)

Updated 11/2/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-21584](http://www.novusbio.com/reviews/destination/NBP2-21584)



**NBP2-21584**

## OCT1 Antibody


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS, 1% BSA, 20% Glycerol
Target Molecular Weight	76 kDa

Product Description	
Description	Novus Biologicals Rabbit OCT1 Antibody (NBP2-21584) is a polyclonal antibody validated for use in IHC, WB, ICC/IF and IP. Anti-OCT1 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5451
Gene Symbol	POU2F1
Species	Human, Rat
Reactivity Notes	Xenopus laevis (81%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human OCT1. The exact sequence is proprietary.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunoprecipitation 1:100-1:500, Immunohistochemistry-Paraffin 1:100-1:1000

**Images**

Western Blot: OCT1 Antibody [NBP2-21584] - Various whole cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with OCT1 antibody diluted at 1:1000.




Immunohistochemistry-Paraffin: OCT1 Antibody [NBP2-21584] - Rat brain. OCT1 antibody diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min.



Immunohistochemistry-Paraffin: OCT1 Antibody [NBP2-21584] - HeLa xenograft, using OCT1 antibody at 1:750 dilution. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



Western Blot: OCT1 Antibody [NBP2-21584] - Wild-type (WT) and OCT1 knockout (KO) HeLa cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with OCT1 antibody diluted at 1:2000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.



OCT1 antibody detects OCT1 protein at nucleus by immunofluorescent analysis.


Sample: HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min.

Green: OCT1 protein stained by OCT1 antibody (NBP2-21584) diluted at 1:500.


Red: phalloidin, a cytoskeleton marker, stained by ( ) diluted at 1:200.

Blue: Hoechst 33342 staining.

Scale bar = 10 um.



OCT1 antibody immunoprecipitates OCT1 protein in IP experiments. IP Sample: 293T whole cell lysate/extract A : Control with 2.5 ug of pre-immune rabbit IgG B : Immunoprecipitation of OCT1 by 2.5 ug of OCT1 antibody (NBP2-21584) 7.5% SDS-PAGE The immunoprecipitated OCT1 protein was detected by OCT1 antibody (NBP2-21584) diluted at 1 : 1000. EasyBlot anti-rabbit IgG (HRP) was used as a secondary reagent.



## Publications

Perovanovic J, Wu Y, Abewe H et al. Oct1 cooperates with the Smad family of transcription factors to promote mesodermal lineage specification *Science signaling* 2023-04-18 [PMID: 37071732]



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-21584**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-21584](http://www.novusbio.com/reviews/submit/NBP2-21584)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

