

Product Datasheet

MAVS Antibody (58N3E1) NBP2-22378

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-22378

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-22378



NBP2-22378

MAVS Antibody (58N3E1)

Product Information

Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	58N3E1
Preservative	0.05% Sodium Azide
Isotype	IgG2b Kappa
Purity	Protein G purified
Buffer	PBS and 0.05% BSA

Product Description

Description	Novus Biologicals Mouse MAVS Antibody (58N3E1) (NBP2-22378) is a monoclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	57506
Gene Symbol	MAVS
Species	Human
Immunogen	Within the range of amino acids 93-298 of human MAVS (isoform 1) were used as the immunogen for this antibody.

Product Application Details

Applications	Western Blot
Recommended Dilutions	Western Blot 0.1-2 ug/ml




Images

Western Blot: MAVS Antibody (58N3E1) [NBP2-22378] - Western blot analysis of HEK293 lysate (35 ug per lane, RIPA buffer) using MAVS antibody (NBP2-22378) at 0.3ug/ml. Band observed at ~65kDa. (Expected MW of 56.5kDa according to NP_065797.2).

250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

Western Blot: MAVS Antibody (58N3E1) [NBP2-22378] - Analysis of A) partial recombinant protein and B) human HeLa lysate using IPS-1 antibody at 0.1 ug/ml and 2 ug/ml, respectively. Goat anti-mouse Ig HRP secondary antibody and PicoTect ECL substrate solution were used for this test.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-22378

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)
NBP2-22378B	MAVS Antibody (58N3E1) [Biotin]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-22378

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

