

Product Datasheet

JunD Antibody - BSA Free

NBP2-22439

Unit Size: 200 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-22439

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-22439



NBP2-22439

JunD Antibody - BSA Free


Product Information	
Unit Size	200 ug
Concentration	LYOPH mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Reconstitution Instructions	Reconstitute with 0.2 ml PBS to a final concentration of 1 mg/ml.
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS

Product Description	
Description	Novus Biologicals Rabbit JunD Antibody - BSA Free (NBP2-22439) is a polyclonal antibody validated for use in WB, Flow, ICC/IF and ChIP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3727
Gene Symbol	JUND
Species	Mouse
Immunogen	Synthetic peptide corresponding to residues G(329) C Q L L P Q H Q V P A Y (341) of mouse JunD.

Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Chromatin Immunoprecipitation (ChIP)
Recommended Dilutions	Western Blot 10 ug/ml, Flow Cytometry 3 - 5 ug / 10 ⁶ cells, Immunocytochemistry/ Immunofluorescence 2 - 3 ug/ml, Chromatin Immunoprecipitation (ChIP) 1-3 ul
Application Notes	WB: Detects an approx. 37 - 41 kDa protein representing JunD from both un-induced and induced murine NIH 3T3 cells.

Images

Chromatin Immunoprecipitation: JunD Antibody [NBP2-22439] - JunD analysis using cross-linked chromatin from 1x10⁶ HTC-IR rat hepatoma cells treated w/ insulin @ 0, 10, & 30 mins. IP using multiplex Matrix ChIP assay, 1 ul/100ul well volume of JunD pAb. Chromatin aliquots from ~1x10⁵ cells used/ChIP pull-down. Quantitative PCR data: 4x using 1ul eluted DNA in 2ul SYBR real-time PCR reactions containing primers to amplify exon-1 - exon-28 of LAMC1. PCR calibration curves made per primer pair from dilution series of sheared total genomic DNA. ChIP quantitation: signal relative to total amount of input chromatin. Results: mean +/- SEM for 3 experiments. Schematic representation of rat LAMC1 locus above data where boxes represent exons (black = translated regions, white = untranslated regions). Lines: zigzag represents an intron, & the straight represents upstream sequence. Regions amplified by LAMC1 primers represented by black bars. Data courtesy of the Innovators Program.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-22439

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-22439

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

